



Sacace
BIOTECHNOLOGIES

CATALOGUE 2016

Molecular diagnostics
kits and instruments

Company portrait

Sacace Biotechnologies s.r.l., is an innovative Italian company, located in **Como (Italy) founded in 2001**. From the very beginning this company was focused on innovation and on establishing an international presence. This multinational presence reinforces our ability to offer our healthcare solutions and to anticipate needs in all regions of the world. The company supplies its products to the international diagnostics and pharmaceutical industry, as well as to hospitals and laboratories through a global network of distributors located in more than **60 countries** all over the world.

Sacace prime objective is to identify and meet the customers' needs. This implies solving their problems and anticipating their future needs by maintaining close contact with them and listening to what they say.

Sacace incorporates all functions required – research & development, manufacturing, logistics, technical support, marketing and sales of an extensive line of molecular biology diagnostic tests that accurately screen for the presence of disease in human and food fields – to provide adequate support to its customers.

The kits produced by Sacace include all the required reagents for nucleic acid extraction, reverse transcription, amplification of specific genomic regions and detection by agarose gel electrophoresis, **Real Time PCR**.

The products are designed to provide ease-of-use, they are **CE marked** in compliance with **Directive 98/79EC** and incorporate the highest quality reagents to ensure consistency, reliability and long shelf life.

We have a dedicated team of experienced research scientists who are committed to researching and developing new products and improved test methods to aid patient diagnosis.

Customer consultation is a key element in our strategy for the development and production of kits that meet market needs and exceed expectations.

Quality assurance is the cornerstone of the company's success. The quality management has obtained **ISO 9001:2008** and **ISO 13485:2012** certifications and all aspects of product design and manufacture are carried out in accordance with these standards.



9001 : 2008



13485 : 2012



CERTIFICATO DEL SISTEMA DI GESTIONE PER LA QUALITÀ QUALITY MANAGEMENT SYSTEM CERTIFICATE

Si dichiara che il sistema di gestione per la Qualità dell'Organizzazione:
We certify that the Quality Management System of the Organization:

Reg. No: 8315 - A

Indirizzo/Address:

Via Scalabrini 44
22100 Como Italia

**SACACE
BIOTECNOLOGIES S.r.l.**

È conforme alla norma/Is in compliance with the standard:

UNI EN ISO 9001:2008
ISO 9001:2008

Per i seguenti prodotti-servizi/For the following products-services:

Sviluppo e produzione di kit diagnostici e per la ricerca scientifica
Development and production of diagnostic and scientific research kits

EA: 12

Il mantenimento della certificazione è soggetto a sorveglianza annuale e subordinato al rispetto dei requisiti essenziali CERMET.
Maintenance of the certification is subject to annual survey and dependent upon the observance of CERMET basic requirements.

La data di rilascio di questo certificato corrisponde alla data di primo rilascio da parte di altro Ente accreditato
The date of issuance of this certificate is the date of first issue by another accredited body

Rilascio certificato/Certificate issuance: 2007-06-22
Ultima modifica/Last modification: 2013-03-08
Prossimo rinnovo/Following renewal: 2016-03-09

Direttore Commerciale e Operativo
Sales and Operations Manager
Giampiero Belcredi

Direttore Generale
General Manager
Rodolfo Trippodo



SQO N° 067A
SQA N° 0190
PRQ N° 069B

FSM N° 004
SSI N° 0090





CERTIFICATO DEL SISTEMA DI GESTIONE PER LA QUALITÀ QUALITY MANAGEMENT SYSTEM CERTIFICATE

Si dichiara che il sistema di gestione per la Qualità dell'Organizzazione:
We certify that the Quality Management System of the Organization:

Reg. No: 8315 – M

Indirizzo/Address:

Via Scalabrini 44
22100 Como Italia

**SACACE
BIOTECNOLOGIES S.r.l.**

È conforme alla norma/Is in compliance with the standard:

UNI CEI EN ISO 13485:2012
ISO 13485:2003

Per i seguenti prodotti-servizi/For the following products-services:

Sviluppo e produzione di kit diagnostici
Development and production of diagnostic kits

EA: 12

Il mantenimento della certificazione è soggetto a sorveglianza annuale e subordinato al rispetto dei requisiti essenziali CERMET.

Maintenance of the certification is subject to annual survey and dependent upon the observance of CERMET basic requirements.

Riferirsi al manuale qualità per i dettagli delle esclusioni ai requisiti della norma UNI CEI EN ISO 13485:2012

Refer to quality manual for details of exclusion of ISO 13485:2003 requirements

La presente certificazione è stata rilasciata in conformità al Regolamento Tecnico Accredia / Sincert RT 20.

This certification has been granted in compliance with the Accredia / Sincert Technical Regulation TR 20.

Per informazioni puntuali e aggiornate circa eventuali variazioni intervenute dello stato della certificazione di cui al presente certificato, si prega di contattare il n° telefonico +39 051.459.3.111 o e-mail: infotorino@cermet.it In case of punctual and updated information about any changes to the certification status, please contact phone number +39 051.459.3.111 or e-mail: infotorino@cermet.it

Rilascio certificato/Certificate issuance: 2010-03-10

Ultima modifica/Last modification: 2013-03-08

Prossimo rinnovo/Following renewal: 2016-03-09

Direttore Commerciale e Operativo
Sales and Operations Manager
Giampiero Belcredi

Direttore Generale
General Manager
Rodolfo Trippodo



SGQ N° 007A
SGA N° 0100
PRD N° 0690

FSM N° 0046
SII N° 0060



Molecular Diagnostics

The analysis of **DNA**, **RNA** and proteins at the molecular level performed in clinical laboratories, known as molecular diagnostics, is our core business.

Modern analyses based on the detection of nucleic acids offer considerable advantages over traditional methods of pathogen detection in humans. These procedures detect viruses, bacteria, and parasites more rapidly and with far greater sensitivity and specificity. At genetic and protein level, the cause of a disease can now be found more precisely, enabling the most suitable therapy to be developed. Likewise, the analysis of an individual's genetic makeup enables physicians to predict the course of certain diseases. Therewith, molecular diagnostics provide modern medicine with the necessary tools for developing completely new, personalized strategies in the battle against many diseases. Hospitals and diagnostics laboratories using these techniques have different needs than customers in other markets. They ask for products which guarantee the highest levels of reliability and maximum speed. The advent of molecular biology diagnostics has particularly revolutionized the diagnosis and treatment of diseases.

The company offers a wide range of assays for **real-time polymerase chain reaction (PCR)** applications and these products are sold either direct to an end user or through a distribution relationship.

These products are 'platform independent' and are used by customers on different platforms like Rotor-Gene™ (Qiagen), LineGeneK™ (Bioer Technologies), iQ5, CFX-96™ (BioRad), SmartCycler™ (Cepheid), Applied Biosystems 7300/7500™, EcoqPCR™ (Illumina), MX3005P™ (Agilent Technologies), SaCycler-96™ (Sacace Biotechnologies).

Examples of the diseases tested for include: **HCV, HBV, HIV, sexually transmitted infections, cardiovascular diseases, the herpetic family of viruses (CMV, EBV, and HSV)** as well as seasonal infectious diseases such as **enterovirus** and **influenza**.

Sacace offers different kits for HCV qualitative, quantitative and genotyping tests.

The company's portfolio consists also of multiplexing assays which allow for the testing of several different pathogens in one single run: for example we have a kit for simultaneous multiplex detection of HCV/HBV/HIV, CMV/EBV/HHV6, Chlamydia/Neisseria/Trichomonas/Mycoplasma and others.

One of the most important products lines of Sacace is sexually transmitted diseases (STD) kits. STD refer to a variety of bacterial, viral and parasitic infections that are acquired through sexual activity.

In light of the epidemic spread of animal diseases such as avian flu or swine flu and the recent international food safety scandals, molecular testing is continuously gaining importance in veterinary medicine and agricultural industry. Our molecular tests enable for reliable and rapid detection of infectious diseases and food-borne pathogens.

Sacace offers also a range of nucleic acid purification products that are used to investigate bacterial and viral infections, furthermore proposes a wide range of applications for bacterial drug resistance detection and **molecular genetics**.



Sacace Biotechnologies Srl


via Scalabrini, 44. – 22100 –Como – Italy Tel +390314892927 Fax +390314892926 VAT: 01294510621
mail: info@sacace.com web: www.sacace.com

Table of Contents

Hepatitis Viruses	9
Hepatitis A (HAV)	9
Hepatitis B (HBV)	10
Hepatitis B Genotyping	11
Hepatitis C (HCV)	12
HCV associated infection kits (IL28)	13
Hepatitis C Virus Genotyping	14
Hepatitis D (HDV)	15
Hepatitis G (HGV)	15
Human Immunodeficiency Virus	16
HIV	16
HIV DNA	16
HIV associated infections	17
HCV/HBV/HIV	17
Human Papilloma Virus	18
Kits for screening of high carcinogenic risk HPV genotypes	18
Kits for typing of high carcinogenic risk HPV genotypes	18
Low carcinogenic risk HPV genotypes	19
Sexually transmitted diseases	20
Chlamydia trachomatis	20
Ureaplasma species and Mycoplasma	20
Neisseria gonorrhoeae	21
Trichomonas vaginalis	22
Sexually transmitted Diseases: Multiplex PCR Kits	22
Candida albicans	23
Treponema pallidum	23
Bacterial vaginosis	24
Herpes Virus infections	25
Cytomegalovirus (CMV)	25
Herpes simplex	25
Epstein-Barr virus (EBV)	26
Human herpesvirus 6 (HHV-6)	26
TORCH Infections	27
Rubella	27
Toxoplasmosis	27
Parvovirus B19 (B19V)	27
Respiratory Infections	28
Influenza	28
Chlamydia & Mycoplasma pneumoniae	29
Pseudomonas aeruginosa	29
Legionella pneumophila	29
Human respiratory syncytial virus (hRSV)	30
Acute Respiratory Viral Infections	30
Coronavirus (MERS+SARS)	30
Bordetella Pertussis	31
Tuberculosis detection and differentiation	31
Neurological infections	32
Enteroviruses	32
Meningitis	33
Intestinal and parasitic infections	34
Acute Intestinal Infections (A.I.I.)	34
Rotavirus/Norovirus/Astrovirus	34
Shigella/Salmonella/Campilobacter	34
Aerobic complex	34
Food Pathogens	35
Helicobacter Pylori	36
Tropical parasitic infections	36

Dangerous infections	37
Lyme disease	37
Brucellosis	37
Dengue virus	37
Pestis	37
Bacillus Anthracis	38
Vibrio Cholerae	38
Leptospira	38
Congo Crimea	38
Tick-Borne diseases	38
West Nile virus	39
Coxiella burnetii	39
Bacterial pyogenic: Streptococcus agalactiae, Streptococcus pyogenes, Pseudomonas aeruginosa	39
Ebola	39
Bacterial Drug resistance	40
Molecular Genetics	41
Coagulation / Fibrinolysis.....	41
Pharmacogenetics	44
Diabetes and Obesity	46
Hepatitis C treatment prognosis	47
Hypertension	48
Myocardial infarction / ischemic stroke	49
Woman's Health / Miscarriage	50
Immunity / Lipid Metabolism	52
Men's Health / HLA	53
Oncological diseases	54
Cancer Therapy / Periodontitis	55
DNA/RNA Purification	56
Instruments	58
SaCycler-96	58
SaMag 12/24	59

Kits Table Legend:



Hepatitis Viruses

Hepatitis Viruses are a group of infectious liver diseases caused by hepatotropic viruses belonging to different families. There are 5 major viruses that cause hepatitis. They make up two groups of hepatitis: enteric (HAV and HEV) and parenteral (HBV, HCV and HDV).

	HEPATITIS A	HEPATITIS B	HEPATITIS C	HEPATITIS D	HEPATITIS E
AGENT	Hepatitis A virus (HAV); single stranded RNA; no envelope	Hepatitis B virus (HBV); double stranded DNA; envelope	Hepatitis C virus (HCV); single stranded RNA; envelope	Hepatitis D virus (HDV); single stranded RNA; envelope from HBV	Hepatitis E virus (HEV); single stranded RNA; no envelope
FAMILY	Picornaviridae	Hepadnaviridae	Flaviviridae	Deltavirus (viroid)	Calicivirus
GENOME	RNA: 7500 nc	DNA: 3200 nc	RNA: 9500 nc	RNA: 1700 nc	RNA: 7500 nc
TRANSMISSION WAY	Fecal-Oral	Parenteral	Parenteral	Parenteral	Fecal-Oral
CHRONIZATION	No	Yes	Yes	Yes	No
INCUBATION TIME	15-50 days	45-160 days	14-180 days	Uncertain	15-50 days
MANIFESTATION OR SYMPTOMS	Mostly subclinical; severe cases: fever, headache, malaise, jaundice	Mostly subclinical; similar to HAV, but fever, headache absent, and often progress to severe liver damage	Similar to HBV	Severe liver damage, high mortality rate	Similar to HAV, but pregnant women may have high mortality rate
VACCINES	Yes	Yes	None	HBV vaccine is protective because coinfection required	None

Hepatitis A (HAV)

The **hepatitis A (HAV)** virus is the enteric infection most widely spread in the world. This is the acute infectious disease of liver transmitted by fecal-oral way, the causative agent of which is hepatitis A virus (HAV) belonging to the family Picornaviridae. Virus hepatitis A is one of the five most economically significant infectious diseases and one of the priority problems of the public healthcare.

Detection of the causative agent RNA by PCR method has significant advantages as related to ELISA and biochemical tests at detection of the virus in blood of contact persons as RNA of the hepatitis A virus manifests itself in the blood on the third week from the moment of contamination and is detected at the average within 20 days after appearance of the disease symptoms. Thus, RNA is the first diagnostic marker detected in the patient blood, occurs earlier than HAV IgM and gives no false negative reactions.

Hepatitis A Virus Kits

TV4-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	HAV Real-TM Qual Complete Real Time Test with Ribo-Sorb extraction kit	R	€	50	5 x10 ² copies/ml
TV4-50FRT C SA, RG, IQ, SC, MX, A, IL, B, LC	HAV Real-TM Qual Complete Real Time Test with Ribo-Virus column extraction kit	R	€	50	5 x10 ² copies/ml
V4-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	HAV Real-TM Qual Real Time Amplification kit with the RNA extraction controls	R	€	50	5 x10 ² copies/ml
V-4-50R	HAV 430	A		55	1 x10 ³ copies/ml

Hepatitis B (HBV)

Hepatitis B virus (HBV) is a widely spread human infection caused by DNA-containing virus of hepatitis B belonging to the family Hepadnaviridae. Transmission of hepatitis B virus results from exposure to infectious blood or body fluids containing blood. Possible forms of transmission include unprotected sexual contact, blood transfusions, re-use of contaminated needles & syringes, and vertical transmission from mother to child during childbirth. The viral hepatitis B presents a serious problem for public healthcare due to its universal spread. At present in accordance with the WHO data the population infected with hepatitis B virus makes 500 million people.

Detection of HBV DNA is used for:

- Early diagnostics of acute viral hepatitis B;
- Detection of latent forms of viral hepatitis B;
- Detection of mutant strains of hepatitis B virus by HBsAg;
- Establishment of diagnosis of chronic viral hepatitis B;
- Monitoring of effectiveness of the antiviral therapy;

ADVANTAGES OF SACACE™ HBV REAL-TM QUANT DX CE-IVD MARKED KIT

Key Features

- Primers and probes in highly conserved 5'-gene (coding the surface antigen of the hepatitis B virus HBsAg) region of the HBV genome.
- **Reagents lyophilized** and aliquoted
- Excellent **sensitivity** of **7 IU/ml** (1 ml input)
- Results expressed directly in International Units
- Use of exogenous **internal control** to check extraction and amplification
- 1 Low and 1 High concentrated extraction positive controls, to be extracted as the samples
- Reagents can be **stored at 4°C** and **shipped at room temperature**
- Very **long shelf life** (1 year)

Advantages


- No need of PCR mix setup
- No possibility of mistake during reagent dispensing
- No problem with storage
- No need for refrigerated transportation
- No possibility of reagent components contamination

Very Easy To Use

No more need to prepare PCR mastermix! The HBV Real-TM Quant DX kit contains **96 ready to use 0.2 ml PCR tubes where you just need to add your extracted viral DNA. That's it!**


The reaction tubes are ready to be transferred into the Real Time PCR thermal cycler.

Excellent sensitivity (LOD)



The Limit Of Detection of 7 IU/ml was determined by testing dilutions of the 3rd WHO International Standard for Hepatitis B Virus for Nucleic Acid Amplification Techniques (NIBSC code: 10/264) prepared in HBV negative human plasma. The results were determined by Probit analysis.

Optimal specificity



Most relevant HBV genotypes (A-B-C-D-H) were tested analyzing dilution series and results showed good and clear sigmoid-shaped fluorescence curves and a Cycle threshold (Ct) less than 35.

Hepatitis B Virus Kits

V5-96/3FRT SA, RG*	HBV Real-TM Quant DX Real Time Amplification kit with positive controls and standards, 96 ready to use lyophilized PCR tubes * validated on SA and RG, but optimized also on iQ,SC,MX,A,B,LC	R	CE	96	Linearity: 7 - 10 ⁸ IU/ml
TV5-100/2FRT SA, RG, iQ, SC,MX, A,IL,B	HBV Real-TM Quant Complete Real Time Test with Ribo-Sorb extraction kit (25 µl Reaction Mix)	R		100	Linearity: 3 x10 ² -10 ⁸ copies/mL
TV5-100/2FRT C SA, RG, iQ, SC,MX, A,IL,B	HBV Real-TM Quant Complete Real Time Test with Ribo-Virus column extraction kit (25 µl Reaction Mix)	R		100	Linearity: 2 x10 ² -10 ⁸ copies/mL
V5-100/2FRT SA, RG, iQ, SC,MX, A,IL,B	HBV Real-TM Quant Real Time Amplification kit with the DNA extraction controls (25 µl Reaction Mix)	R		100	Linearity: 50 -10 ⁸ copies/mL
TV5-100FRT SA, RG, iQ, SC,MX, A,IL,B,LC	HBV Real-TM Qual Complete Real Time Test with Ribo-Sorb extraction kit (25 µl Reaction Mix)	R		100	2 x10 ² copies/ml
TV5-100FRT C SA, RG, iQ, SC,MX, A,IL,B,LC	HBV Real-TM Qual Complete Real Time Test with Ribo-Virus column extraction kit (25 µl Reaction Mix)	R		100	1 x10 ² copies/ml
V5-100FRT SA, RG, iQ, SC,MX, A,IL,B,LC	HBV Real-TM Qual Real Time Amplification kit with the DNA extraction controls (25 µl Reaction Mix)	R		100	50 copies/ml
TV5-64FRT SA, RG, iQ, SC,MX, A,IL,B,LC	HBV Real-TM Qual Complete Real Time Test with Ribo-Sorb extraction kit (50 µl Reaction Mix)	R		64	1 x10 ² copies/ml
V-5-100R	HBV 470/770 IC	A		110	5 x10 ² copies/ml

Hepatitis B Genotyping

Hepatitis B virus (HBV) infects nearly two billion people worldwide. The hepatitis B virus (HBV) is currently categorized into eight genotypes (A to H). Genotypes have been found to be geographically distributed. Numerous studies have investigated the clinical implications of HBV genotypes to disease severity, response to IFN, disease chronicity and hepatocellular carcinoma (HCC).

Hepatitis B Genotyping Kits

R5-Gen SA, RG, iQ, SC,MX, A,IL,B,LC	HBV Genotype A, B, C, D Real-TM Real Time Amplification kit	R		50	
--	---	---	--	----	--

Hepatitis C (HCV)

Hepatitis C virus (HCV) is RNA-containing, hepatotropic virus belonging to the Flaviviridae family. Contamination with hepatitis C virus occurs at direct entering of the virus in blood (at parenteral interventions or during blood transfusions). Most people with acute HCV infection are asymptomatic or have mild symptoms (fatigue, nausea, jaundice) but they are unable to clear the virus and in approximately 80% of cases this leads to chronic infection. In 15 to 20% of patients chronic HCV infection progresses at a variable rate to cirrhosis, with a 1 to 4% annual risk of developing hepatocellular carcinoma.

ADVANTAGES OF SACACE™ HCV REAL-TM QUANT DX CE-IVD MARKED KIT

Key Features

- Primers and probes in highly conserved 5' UTR region of the HCV genome
- **Reagents lyophilized** and aliquoted
- Excellent **sensitivity** of **13 IU/ml** (1 ml input)
- Results expressed directly in International Units
- Use of exogenous **internal control** (RNA) to check extraction and amplification
- 1 Low and 1 High concentrated extraction positive controls, to be extracted as the samples
- Reagents can be **stored at 4°C** and **shipped at room temperature**
- Very **long shelf life** (1 year)

Advantages


- No need of PCR mix setup
- No possibility of mistake during reagent dispensing
- No problem with storage
- No need for refrigerated transportation
- No possibility of reagent components contamination

Very Easy To Use

No more need to prepare PCR mastermix! The HCV Real-TM Quant DX kit contains **96 ready to use 0.2 ml PCR tubes where you just need to add your extracted viral RNA. That's it!**

The reaction tubes are ready to be transferred into the Real Time PCR thermal cycler.


Optimal specificity



Most relevant HCV genotypes (1-2-3-4-5) were tested analyzing dilution series and results showed good and clear sigmoid-shaped fluorescence curves and a Cycle threshold (Ct) less than 35.



Wide linear range



The linear range of the HCV Real-TM Quant DXkit has been determined by analyzing a dilution series (8,00 log IU/ml to 1,00 log IU/ml) of an HCV synthetic quantitative standard calibrated against the 4th WHO International HCV RNA Standard.

Hepatitis C Virus Kits

V1-96/3FRT SA, RG*	HCV Real-TM Quant DX Real Time Amplification kit with positive controls and standards, 96 ready to use lyophilized PCR tubes * validated on SA and RG, but optimized also on iQ,SC,MX,A,B,LC	R	CE	96	Linearity: 13 - 10 ⁸ IU/ml
TV1-100/2FRT SA, RG, iQ, SC,MX, A,IL,B,	HCV Real-TM Quant Complete Real Time Test with Ribo-Sorb extraction kit (25 µl Reaction Mix)	R		100	Linearity: 3 x10 ² - 5 x10 ⁷ IU/mL
V1-100/2FRT SA, RG, iQ, SC,MX, A,IL,B	HCV Real-TM Quant Real Time PCR kit with the RNA extraction controls (25 µl Reaction Mix)	R		100	Linearity: 50 - 5 x10 ⁷ IU/mL
TV1-100/2FRT C SA, RG, iQ, SC,MX, A,IL,B	HCV Real-TM Quant Real Time PCR kit with Ribo-Virus column extraction kit (25 µl Reaction Mix)	R		100	Linearity: 2 x10 ² - 5 x10 ⁷ IU/mL
TV1-100FRT SA, RG, iQ, SC,MX, A,IL,B,LC	HCV Real-TM Qual Complete Real Time Test with Ribo-Sorb extraction kit (25 µl Reaction Mix)	R		100	2 x10 ² IU/mL
TV1-100FRT C SA, RG, iQ, SC,MX, A,IL,B,LC	HCV Real-TM Qual Complete Real Time Test with Ribo-Virus column kit (25 µl Reaction Mix)	R		100	1 x10 ² IU/mL
V1-100FRT SA, RG, iQ, SC,MX, A,IL,B,LC	HCV Real-TM Qual Real Time PCR kit with the RNA extraction controls (25 µl Reaction Mix)	R		100	50 IU/mL
V-1-100R	HCV - 240/440 IC (RT-PCR)	A		110	5 x10 ² IU/mL

HCV associated infection kits

Interleukin-28 (IL28) is a cytokine that plays a role in immune defense against viruses. IL28B belongs to the type III interferon family of cytokines. Its classification as interferon is due to its ability to induce an antiviral state. Polymorphisms in the IL28B gene region are important in predicting outcome following therapy for chronic hepatitis C virus (HCV) infection.

Combined therapy IFN pegylated (PEG-IFN) and ribavirin (RBV) is the current standard therapy against HCV infection and to know in detail the polymorphism in IL28B gene region of patients infected with HCV can be an important component of the decision to initiate treatment with PEG-IFN and RBV.

HCV associated infection kits


R05-100FRT SA, RG, iQ, MX, A, B	IL28B rs17 / rs60 Real-TM Real Time amplification kit	R		100	
------------------------------------	---	---	--	-----	--

Hepatitis C Virus Genotyping

HCV is classified into eleven major genotypes (designated 1-11), many subtypes (designated a, b, c, etc.), and about 100 different strains (numbered 1,2,3, etc.) based on the genomic sequence heterogeneity. Genotypes 1-3 have a worldwide distribution. Types 1a and 1b are the most common, accounting for about 60% of global infections. They predominate in Northern Europe and North America, and in Southern and Eastern Europe and Japan, respectively. Type 2 is less frequently represented than type 1. Type 3 is endemic in south-east Asia and is variably distributed in different countries. Genotype 4 is principally found in the Middle East, Egypt, and central Africa. The determination of the infecting genotype is important for the prediction of response to antiviral treatment: genotype 1 and 4 are generally associated with a poor response to interferon alone, whereas genotypes 2 and 3 are associated with more favourable responses. At patients with subtype 1b the disease progresses to a chronic condition 90 % of cases, in that time as with genotypes 2 and 3b in 33-50 %. In a number of works it is mentioned, that infection with 1b genotype have heavier current of disease with development of a cirrhosis and hepatocarcinoma.

The International Consensus European Association for the Study of the Liver (EASL) recommends before beginning of antiviral therapies to carry out a liver biopsies and to determine HCV genotype.

Global distribution of HCV genotypes



Hepatitis C Virus Genotyping Kits

R1-Gen-4X SA, RG, IQ, MX, SC, A, B	HCV 1/2/3 Genotype Real-TM Real Time PCR kit with the RNA extraction controls	R	48	1 x10 ² IU/mL
R1-Gen-6 SA, RG, IQ, SC, MX, A, IL, B	HCV Genotype Plus (1a,1b, 2, 3a, 4, 5a, 6) Real-TM Real Time PCR kit with the RNA extraction controls	R	48	1 x10 ² IU/mL
V-1-G/50R	HCV - Genotype 1a, 1b, 2, 3 (RT-PCR)	A	55	1 x10 ³ IU/mL

Hepatitis D (HDV)

Hepatitis D virus (HDV) is RNA containing, hepatotropic viroid (uncompleted virus) belonging to Deltavirus family. HDV needs helper function of hepatitis B virus that provides to HDV proteins of the superficial membrane (HBsAg) that's why HDV can replicate itself only in presence of HBV. Transmission of HDV can occur either via simultaneous infection with HBV (coinfection) or via infection of an individual previously infected with HBV (superinfection). Both superinfection and coinfection with HDV results in more severe complications compared to infection with HBV alone. These complications include a greater likelihood of experiencing liver failure in acute infections and a rapid progression to liver cirrhosis, with an increased chance of developing liver cancer in chronic infections. In combination with hepatitis B virus, hepatitis D has the highest mortality rate of all the hepatitis infections of 20%.

Detection of HDV RNA by PCR allows detection of the causative agent in the period of introduction of infection before seroconversion, which is very important for early diagnostics.

Hepatitis D Virus Kits

TV3-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	HDV Real-TM Qual Complete Real Time Test with Ribo-Sorb extraction kit	R	50	5 x10 ² copies/mL
TV3-50FRT C SA, RG, IQ, SC, MX, A, IL, B, LC	HDV Real-TM Qual Complete Real Time Test with Ribo-Virus column kit	R	50	3 x10 ² copies/mL
V3-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	HDV Real-TM Qual Real Time PCR kit with the DNA extraction controls	R	50	1 x10 ² copies/mL
V3-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	HDV Real-TM Qual Real Time PCR kit with the DNA extraction controls	R	100	1 x10 ² copies/mL
V3-100/2FRT SA, RG, IQ, SC, MX, A, B	HDV Real-TM Quant Real Time PCR kit with the DNA extraction controls	R	100	1 x10 ² copies/mL
TV56-100FRT SA, RG, IQ, SC, MX, A, B	HBV/HDV Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R	100	2 x10 ² /5-10 ² copies/mL
TV56-100FRT C SA, RG, IQ, SC, MX, A, B	HBV/HDV Real-TM Complete Real Time Test with Ribo-Virus column kit	R	100	1 x10 ² /3-10 ² copies/mL
V56-100FRT SA, RG, IQ, SC, MX, A, B	HBV/HDV Real-TM Real Time PCR kit with the DNA/RNA extraction controls	R	100	20/100 copies/mL
V-3-100R	HDV 255/500 IC (RT-PCR)	A	110	2 x10 ³ copies/mL

Hepatitis G (HGV)

Hepatitis G virus (HGV) is another virus causing post-transfusion hepatitis. The same as hepatitis C virus, hepatitis G virus belongs to the flaviviruses. The hepatitis G virus is detected with the help of PCR (serological methods are less reliable). It's detected in 1.5 percent of donors and in some patients with acute fulminant and chronic hepatitis. The coinfection of hepatitis G with hepatitis B, C or D is detected frequently.

Hepatitis G Virus Kits

TV2-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	HGV Real-TM Real Time PCR kit with Ribo-Sorb extraction kit	R	50	5 x10 ² copies/mL
V2-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	HGV Real-TM Real Time PCR kit with the DNA extraction controls	R	50	3 x10 ² copies/mL
V-2-100R	HGV 340/625 IC (RT-PCR)	A	110	2 x10 ³ copies/mL

Human immunodeficiency virus (HIV)

HIV is a lentivirus (a member of the retrovirus family) differentiated on structural and antigenic properties into two virus types: HIV-1 and HIV-2. HIV-2 occurs considerably less often than HIV-1. In accordance with 1991 Nomenclature, there are three independent HIV-1 groups: «M» (main); «O» (outlier); «N» (non-V/non-O). Groups O and N are less widely spread and occur in African countries population. Group M includes 11 subtypes: A1, A2, B, C, D, F1, F2, G, H, J, K.

Transmission ways of the virus are very important for the virus spread. HIV is transmitted by three ways: at heterosexual and homosexual intercourse, parenteral with blood and blood products and vertically: from the infected mother to the child by an intrauterine way, during the child delivery or soon after the childbirth at breast feeding.

This method has a lot of advantages:

- detection of virus DNA/RNA allows reducing the length of the “serological window”;
- PCR is an indispensable approach for HIV-diagnostics in children born from HIV-infected mothers;
- determination of HIV RNA in the blood plasma (viral load) is an obligatory procedure to monitoring of the therapy effectiveness

ADVANTAGES OF SACACE™ HIV REAL-TM QUANT KIT

- Application of primers and probes in the most conservative area of the HIV-1 polymerase gene that allow effective detection of the majority of HIV-1 subtypes.
- Use of the Quantitative Internal Control (concentration reported in Data Card) which allows not only to monitor the extraction procedure and to check possible PCR inhibition but also to verify possible losses of the RNA during extraction procedure thus enabling to calculate precisely the HIV viral load.
- Presence in the reagents supplied with the kit of two positive controls of the extraction: Pos1 – low viral load and Pos2 – medium viral load that are quantitatively described in Data Card and allow quality control of the conducted analysis.
- Use of Quantitative Standards for HIV RNA and HIV IC enabling to calculate precisely the HIV viral load.
- The reagent kit possesses a wider linear range of measurements (from 25 to 5×10^6 copies/ml).

HIV RNA Quant Kits

V0-96/3FRT SA, RG*	HIV Real-TM Quant Dx Real Time PCR Test with positive controls and standards (96 ready to use lyophilized tubes - 50 µl Reaction Mix) <small>* validated on SA and RG, but optimized also on iQ,SC,MX,A,B,LC</small>	R	€€	96	Linearity: 48 - 1×10^7 IU/mL
TR-VM-100FRT SA, RG, iQ, MX, A,IL,B	HIV RNA Real-TM Quant Complete Real Time Test with Ribo-Sorb extraction kit (25 µl Reaction Mix)	R		100	Linearity: $2,5 \times 10^2$ - 5×10^6 copies/mL
TR-VM-100FRT C SA, RG, iQ, MX, A,IL,B	HIV RNA Real-TM Quant Complete Real Time Test with Ribo-Virus column kit (25 µl Reaction Mix)	R		100	Linearity: $2,0 \times 10^2$ - 5×10^6 copies/mL
R-VM-100FRT SA, RG, iQ, MX, A,IL,B	HIV RNA Real-TM Quant Real Time PCR kit with the RNA extraction controls (25 µl Reaction Mix)	R		100	Linearity: 50- 5×10^6 copies/mL

HIV DNA

HIV-infection diagnostics in children born from HIV-infected mothers is difficult due to the fact that mother's antibodies to HIV persist in such children's blood for a long time. The problem of earlier HIV-infection diagnostics in newborns was solved with development of molecular-genetic methods that allow detection of HIV genome fragments in the peripheral blood at early infection stages.

HIV DNA kit

TR-V1-D SA, RG, iQ, MX, SC, A	HIV DNA Real-TM Qual Complete Real Time Test with Hemo-Sorb extraction kit	R		100	1×10^2 copies/ml
----------------------------------	--	---	--	-----	---------------------------

HIV Infection associated Kits

Hypersensitivity reaction to abacavir is strongly associated with the presence of the HLA-B*5701 allele. Abacavir (Ziagen, also in the combination pills Kivexa and Trizivir) is a potent antiretroviral drug that is a popular choice for first-line antiretroviral HIV therapy. Its main disadvantage is a hypersensitivity reaction that occurs in between 5% - 8% of patients treated with this drug. **HLA-B*5701 Real-TM** test can predict who will develop a severe allergic reaction to the **anti-HIV drug abacavir** as the presence of HLA-B*5701 is significantly associated with an abacavir hypersensitivity.

HIV Infection associated kits

H53-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	HLA B*5701 Real-TM Real Time Amplification Kit	R	€	100	1 x10 ³ cells/ml
--	--	---	---	-----	-----------------------------

Pneumocystis pneumonia (PCP) or pneumocystosis is a form of pneumonia, caused by the yeast-like fungus (which had previously been erroneously classified as a protozoan) **Pneumocystis jirovecii (carinii)**. P jiroveci is now one of several organisms known to cause life-threatening opportunistic infections in patients with advanced HIV infection worldwide.

TP2-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Pneumocystis jirovecii (carinii) Real-TM Complete Real Time Test with DNA/RNA Prep extraction kit	R	€	50	5 x10 ² copies/ml
P2-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Pneumocystis jirovecii (carinii) Real-TM Real Time Amplification kit	R	€	50	5 x10 ² copies/ml

Cryptococcosis, caused by **Cryptococcus neoformans**, is the most common fungal disease in HIV infected persons and it is the AIDS-defining illness in 60-70% of HIV infected patients.

F4-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Cryptococcus neoformans Real-TM Real Time Amplification Kit	R	€	100	
---	---	---	---	-----	--

HCV/HBV/HIV Real-TM

Transfusion-associated transmission risk of infectious diseases has been reported worldwide. For screening of blood donations in order to reduce the residual risk of transmission of bloodborne viruses, viral nucleic acid testing (NAT), has been introduced by the European plasma industry in 1995, and subsequently introduced for blood donations in many countries. NAT was implemented to reinforce the safety of the blood supply; it can detect acute viral infections during the 'window period', that are not detected by the serological screening methods. Current NAT procedures usually demand pooling of blood donation samples due to the format of the employed platforms.

ADVANTAGES OF SACACE™ HCV/HBV/HIV REAL-TM KIT

- Simultaneously amplification (multiplex) in 1 PCR tube of nucleic acids from HIV, HCV, HBV;
- Separate real-time detection and differentiation of nucleic acids from HIV, HCV and HBV on different channels (FAM – HCV, JOE/HEX/Cy3 – HIV, ROX – HBV, Cy5 – internal control);
- Optimization on different equipments;
- Possibility of pooling (5-10 samples in pool format is recommended);
- High sensitivity*: **HCV RNA** – 10 IU/ml; **HBV DNA** – 5 IU/ml; **HIV RNA** – 20 copies/ml.

* values obtained using the "Magno-Virus" extraction kit (Sacace REF K-2-16/1000)

Hepatitis C / Hepatitis B / HIV multiplex detection kits

V50-100FRT SA, RG, IQ, MX, A, SC	HCV/HBV/HIV Real-TM Real Time PCR Test with RNA/DNA extraction controls	R		100	10/5/20 IU/mL
V62-100FRT RG, SA	HCV/HBV/HIV1/HIV2 Real-TM Real Time PCR Test	R		100	10/5/20 IU/mL

Human Papilloma Virus

Cervical cancer (CC) is one of the most widely spread oncological pathologies that ranks second by the incidence in women in the world. Each year about 600 thousand of new CC cases are registered in the world with more than 250 thousand lethal outcomes. The virus nature of this cancer is confirmed by the World Health Organization and HPV is detected practically in 100 percent of cases of cervical precancer and cancer. Based on the frequency of detection of HPV genotypes from different grades of Cervical Intraepithelial Neoplasia (CIN Grades I – III), HPV genotypes are subdivided into High-risk HPV types (16, 18, 31 and 45), Intermediate-risk types (33, 35, 39, 51, 52, 56, 58, 59, and 68), and Low-risk types (6, 11, 42-44).

Owing to the fact that the cervical cancer (CC) has a long development period and a fail-safe recognizable pre-clinic phase there's a possibility to detect and prevent the disease on its early stage.

Kits for screening of high carcinogenic risk HPV genotypes

The **NEW** kit **HPV 14 Screening & 16,18,45 Typing Real-TM Quant** is an *in vitro* Real Time amplification test for quantitative detection and genotyping of HPV 16,18,45 and simultaneous quantitative detection of HPV types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68 (total 14 genotypes detected). It is known, that the parameter of viral load has a prognostic value and the viral load less than 10^5 HPV genomic equivalents in the swab or 10^3 genomic equivalents for 10^5 cells is considered as insignificant and indicates the presence of transitory infection, however such level of load may have a value only in cases of treatment monitoring. Viral load of more than 10^5 genomic equivalents for 10^5 cells is considered to be important with high significance and indicates the existence of dysplastic changes or high risk of their occurrence. Quantitative detection of viral load allows to evaluate the character of the infection and to make a forecast concerning the stage of the disease.

HPV High Risk Screen Kits


V31-100/F FRT SA, RG, B, iQ, MX, A	HPV 14 Screening & 16,18,45 Typing Real-TM Quant NEW Real Time Amplification kit	R	€€	100	5×10^2 copies/ml
V31-100/2FRT SA, RG, iQ, MX, A, B	HPV High Risk Screen Real-TM Quant Real Time Amplification kit	R	€€	100	5×10^2 copies/ml
TV31-100/2FRT 2X SA, RG, iQ, SC, MX, A,IL,B	HPV High Risk Screen Real-TM 2X Quant (2 channels) Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	5×10^2 copies/ml
V31-100/2FRT 2X SA, RG, iQ, SC, MX, A,IL,B	HPV High Risk Screen Real-TM 2X Quant (2 channels) Real Time Amplification kit	R	€€	100	5×10^2 copies/ml
V-26-50F	HPV High Risk Screen	A		55	1×10^3 copies/ml
V-26-100F	HPV High Risk Screen	A	€€	110	1×10^3 copies/ml

Kits for typing of high carcinogenic risk HPV genotypes


The **NEW** **HPV Genotypes 14 Real-TM** kit is an *in vitro* multiplex Real Time amplification test for qualitative detection and genotyping of up to 14 genotypes of *Human Papillomavirus* (16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68) in urogenital swabs and biopsies.

The **HPV Genotypes 14 Real-TM** kit is based on two major processes: isolation of DNA from specimens and multiplex Real Time amplification of 4 tubes for each sample. The test uses primers directed against regions of *HPV* types and β -globine gene used as Internal Control. If the swab is not correctly prepared (high quantity of mucous or insufficient quantity of epithelial cells) the Internal Control will not be detected.







FAM/Green channel: HPV genotypes 16, 39, 33 and 58



JOE/HEX channel: HPV genotypes 31, 45, 35 and 52



ROX channel: HPV genotypes 18, 59, 56 and 51



Cy5 channel: HPV genotypes 66, 68 and β -Globine

HPV High Risk Typing Kits

V67-100FRT SA, RG, MX, IQ, A, IL, B	HPV genotypes 14 Real-TM Quant NEW Real Time Amplification kit for detection and quantification of high risk genotypes 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68	R	€	100	5×10^2 copies/ml
TV12-100FRT SA, RG, SC, MX, IQ, A, IL, B	HPV 16/18 Real-TM Quant Complete Real Time Test with DNA-Sorb-A extraction kit	R	€	100	5×10^2 copies/ml
V12-100FRT SA, RG, IQ, SC, MX, A, IL, B	HPV 16/18 Real-TM Quant Real Time Amplification kit	R	€	100	5×10^2 copies/ml
V-25-50F	HPV High Risk Typing	A	€	55	1×10^3 copies/ml
V-12-100R	Human Papillomavirus 16/18	A	€	110	5×10^2 copies/ml
V-13-100R	Human Papillomavirus 31/33	A		110	5×10^2 copies/ml
V-14-100R	Human Papillomavirus 35/45	A		110	5×10^2 copies/ml

Low carcinogenic risk

A group of **low carcinogenic risk HPV** is represented by more than 12 genotypes and they are called “low risk” because they cannot cause cervical cancer. Sometimes low-risk HPV types can cause visible changes in the genital area, called genital warts. Genital warts are growths or bumps in the genital areas of men and women. They usually are painless. They may be raised, flat, small or large, and single or multiple. Among low risk HPV the genotypes **6** and **11** are of greatest importance as they are responsible for the overwhelming amount of low-carcinogenic pointed condylomas of genital organs and for more than 90 percent of cases of condylomatosis of the larynx in children.

HPV Low Risk Typing Kits

TV11-100FRT SA, RG, IQ, SC, MX, A, IL, B	HPV 6/11 Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€	100	5×10^2 copies/ml
V11-100FRT SA, RG, IQ, SC, MX, A, IL, B	HPV 6/11 Real-TM Real Time Amplification kit	R	€	100	5×10^2 copies/ml
V-11-50F	HPV 6/11	A	€	55	5×10^2 copies/ml


Sexually Transmitted Diseases

STDs (sexually transmitted diseases) refer to a variety of bacterial, viral and parasitic infections that are acquired through sexual activity. Common STDs include: chlamydia, gonorrhea, herpes, HIV, HPV, syphilis, gardnerella, mycoplasma and trichomoniasis. Many STDs affecting women show no early signs or symptoms. As a result, they go undetected and untreated until complications arise. The consequences of untreated STDs are often more serious in women, including: infertility, tubal pregnancy, chronic pain, cervical cancer and other complications. Early screening, diagnosis, counseling and treatment can stop the spread of STDs.

Enzyme immunoassay(EIA) is the diagnostic method most commonly used for the laboratory diagnosis of STD infections, but EIA has a lower detection limit of 10000 elementary bodies and thus lacks sensitivity required for a screening assay, especially in asymptomatic men. Culture has been the “gold standard” for the diagnosis of many STD and have high sensitivity and specificity however due to a slow-growing tendency it takes 2-3 days to get a result and also requires an invasively taken specimen. Nucleic acid based amplification assays using polymerase chain reaction (PCR) have a lower detection limit of one to 10 elementary bodies and specificities comparable with culture. They also offer all the advantages of non-culture tests in terms of specimen transport, batching, and rapid processing time of approximately 2-3 hours . The improved sensitivity of these assays allows the use of non-invasive specimens such as first catch urine (FCU) specimens. PCR tests using FCU specimens have been shown to have sensitivities ranging from 87% to 97% for men and 82% to 93% for women with specificities of 98–100%.

Chlamydia trachomatis

C. trachomatis can be differentiated into 18 serovars (serologically variant strains) based on monoclonal antibody–based typing assays. Serovars A, B, Ba, and C are associated with trachoma (a serious eye disease that can lead to blindness), serovars D-K are associated with genital tract infections, and L1-L3 are associated with lymphogranuloma venereum (LGV). *Chlamydia trachomatis* is an obligate intracellular pathogen (i.e. the bacterium lives within human cells) and can cause numerous disease states in both men and women. Both sexes can display urethritis, proctitis (rectal disease and bleeding), trachoma, and infertility.



Chlamydia trachomatis Kits

TB1-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Chlamydia trachomatis Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R		100	5 x10 ² copies/ml
B1-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Chlamydia trachomatis Real-TM Real Time Amplification kit	R	€	100	5 x10 ² copies/ml
B1-100/2FRT SA, RG, IQ, SC, MX, A, IL, B	Chlamydia trachomatis Quant Real-TM Real Time Amplification kit	R		100	2,5 x10 ² copies/ml
B86-100FRT SA, RG, IQ, SC, MX, A	Chlamydia trachomatis A, B, C Typing Real-TM Real Time Amplification kit	R		100	1 x10 ³ copies/ml
B-1-50R	Chlamydia trachomatis 330/740 IC	A		55	1 x10 ³ copies/ml
B-1-100R	Chlamydia trachomatis 330/740 IC	A		110	1 x10 ³ copies/ml

Ureaplasma species and Mycoplasma

Ureaplasma species and Mycoplasma are causes of nonchlamydial nongonococcal urethritis. Mycoplasma species do not cause vaginitis, but they may proliferate in patients with bacterial vaginosis and may contribute to the condition. *Ureaplasma species* can cause placental inflammation and may invade the amniotic sac early, causing persistent infection and adverse pregnancy outcomes, including premature birth.

Ureaplasma Kits

TB2-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Ureaplasma species Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€	100	5 x10 ² copies/ml
B2-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Ureaplasma species Real-TM Real Time Amplification kit	R	€	100	5 x10 ² copies/ml

Ureaplasma Kits

B2-100FRT Q SA, RG, IQ, SC, MX, A, IL, B	Ureaplasma species Quant Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
TB19-100FRT SA, RG, IQ, SC, MX, A, IL, B	Ureaplasma parvum/urealyticum Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	1 x10 ³ copies/ml
B19-100FRT SA, RG, IQ, SC, MX, A, IL, B	Ureaplasma parvum/urealyticum Real-TM Real Time Amplification kit	R	€€	100	1 x10 ³ copies/ml
B19-100FRT Q SA, RG, IQ, SC, MX, A, IL, B	Ureaplasma parvum/urealyticum Quant Real-TM Real Time Amplification kit	R	€€	100	1 x10 ³ copies/ml
B-2-100R	Ureaplasma urealyticum 450/750 IC	A	€€	110	1 x10 ³ copies/ml

Mycoplasma Kits

TB3-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Mycoplasma hominis Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	5 x10 ² copies/ml
B3-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Mycoplasma hominis Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
B3-100FRT Q SA, RG, IQ, SC, MX, A, IL, B	Mycoplasma hominis Quant Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
TB4-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Mycoplasma genitalium Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	5 x10 ² copies/ml
B4-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Mycoplasma genitalium Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
TB4-100FRT Q SA, RG, IQ, SC, MX, A, IL, B	Mycoplasma genitalium Quant Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	5 x10 ² copies/ml
B4-100FRT Q SA, RG, IQ, SC, MX, A, IL, B	Mycoplasma genitalium Quant Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
B-3-100R	Mycoplasma hominis 330/550 IC	A	€€	110	1 x10 ³ copies/ml
B-4-100R	Mycoplasma genitalium 280/550 IC	A	€€	110	1 x10 ³ copies/ml

Neisseria gonorrhoeae

Gonorrhea, which is caused by *Neisseria gonorrhoeae*, is an important public health problem and is the most common reportable infectious disease. An estimated 700,000 new gonococcal infections occur annually in the United States. Gonorrhea is most frequently spread during sexual contact. However, it can also be transmitted from the mother's genital tract to the newborn during birth, causing ophthalmia neonatorum and systemic neonatal infection. The incubation period is usually 2-8 days.

Neisseria gonorrhoeae Kits

TB5-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Neisseria gonorrhoeae Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	5 x10 ² copies/ml
B5-100FRT RG, IQ, SC, MX, A, IL, B, LC	Neisseria gonorrhoeae Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
B204-100FRT RG, IQ, SC, MX, A, IL, B	Neisseria gonorrhoeae Quant Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
B-5-100R	Neisseria gonorrhoeae 370/660 IC	A		80+30	1 x10 ³ copies/ml

Trichomonas vaginalis

Trichomonas vaginalis trophozoite is oval as well as flagellated. Five flagella arise near the cytostome; four of these immediately extend outside the cell together, while the fifth flagellum wraps backwards along the surface of the organism. Trichomoniasis is the most common, curable sexually transmitted disease in the world. It is also one of the three most common vaginal infections in women. Trichomoniasis is caused by a one-celled parasite, Trichomonas vaginalis. Trichomoniasis affects both women and men. The most common location of infection in women is the vagina, and in men it is the urethra. In women, the symptoms of trichomoniasis may include yellow-green vaginal discharge, fishy odor, pain during urination and sexual intercourse, and genital itching or irritation. Men usually do not show trichomoniasis symptoms, but some may experience discharge from the penis or burning during urination or ejaculation.



Trichomonas vaginalis

TB6-100FRT SA, RG, iQ, SC, MX, A, IL, B, LC	Trichomonas vaginalis Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€	100	5 x10 ² copies/ml
B6-100FRT SA, RG, iQ, SC, MX, A, IL, B, LC	Trichomonas vaginalis Real-TM Real Time Amplification kit	R	€	100	5 x10 ² copies/ml
B-6-100R	Trichomonas vaginalis 240/520 IC	A		110	1 x10 ³ copies/ml

Sexually Transmitted Disease: Multiplex PCR Kits

TB43-100FRT SA, RG, iQ, SC, MX, A, B	Chlamydia trachomatis/Ureaplasma/M.hominis Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R		100	5 x10 ² copies/ml
B43-100FRT SA, RG, iQ, MX, SC, A, B	Chlamydia trachomatis/Ureaplasma/M.hominis Real-TM Real Time Amplification kit	R	€	100	5 x10 ² copies/ml
TB46-100FRT SA, RG, iQ, MX, SC, A, B	Chl. trachomatis/Ureaplasma/M.genitalium Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R		100	5 x10 ² copies/ml
B46-100FRT SA, RG, iQ, MX, SC, A, B	Chl. trachomatis/Ureaplasma/M.genitalium Real-TM Real Time Amplification kit	R	€	100	5 x10 ² copies/ml
TB60-100FRT RG, SA-5, B	C.trachomatis/Ureapl./M.hominis/M.genitalium Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R		100	5 x10 ² copies/ml
B60-100FRT RG, SA-5, B	C.trachomatis/Ureapl./M.hominis/M.genitalium Real-TM Real Time Amplification kit	R	€	100	5 x10 ² copies/ml
TB61-100FRT RG, SA-5, B	N.gonor./C.trachomatis/T.vaginalis/M.genitalium Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R		100	5 x10 ² copies/ml
B61-100FRT RG, SA-5, B	N.gonor./C.trachomatis/T.vaginalis/M.genitalium Real-TM Real Time Amplification kit	R	€	100	5 x10 ² copies/ml
TB65-100FRT SA, RG, iQ, SC, MX, A, B	T. vaginalis/N.gonorrhoeae Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€	100	5 x10 ² copies/ml
B65-100FRT SA, RG, iQ, MX, SC, A, B	T. vaginalis/N.gonorrhoeae Real-TM Real Time Amplification kit	R	€	100	5 x10 ² copies/ml
B75-100FRT Q SA, RG, iQ, MX, SC, A, B	Ur. parvum/Ur.urealyticum/M.hominis Quant Real-TM Real Time Amplification kit	R	€	100	5 x10 ² copies/ml
B83-100FRT SA, RG, iQ, MX, SC, A, B	Tr. vaginalis/N.gonorrhoeae/Chl.trachomatis Real-TM Real Time Amplification kit	R		100	5 x10 ² copies/m
TB67-100FRT SA, RG, iQ, MX, SC, A, B	Chl. trachomatis/N.gonorrhoeae/M.genitalium Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R		100	5 x10 ² copies/m
B67-100FRT SA, RG, iQ, MX, SC, A, B	Chl. trachomatis/N.gonorrhoeae/M.genitalium Real-TM Real Time Amplification kit	R	€	100	5 x10 ² copies/m
TB87-100FRT SA, RG, iQ, MX, SC, A, B	HSV I / HSV II / Treponema pallidum Complete Real Time Amplification kit with DNA-Sorb-B extraction kit	R		100	5 x10 ² copies/ml
B87-100FRT SA, RG, iQ, MX, SC, A, B	HSV I / HSV II / Treponema pallidum Real Time Amplification kit	R		100	5 x10 ² copies/ml

Candida albicans


Candida albicans is a diploid fungus (a form of yeast) and a causal agent of opportunistic oral and genital infections in humans. *C. albicans* is commensal and is among the gut flora, the many organisms that live in the human mouth and gastrointestinal tract. Under normal circumstances, *C. albicans* lives in 80% of the human population with no harmful effects, although overgrowth results in candidiasis. Candidiasis is often observed in immunocompromised individuals such as HIV-positive patients but may also occur in the blood and in the genital tract. Candidiasis, also known as “thrush”, is a common condition, usually easily cured in people who are not immunocompromised. To infect host tissue, the usual unicellular yeast-like form of *C. albicans* reacts to environmental cues and switches into an invasive, multicellular filamentous forms.

Candida albicans Kits

TF1-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Candida albicans Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	1 x10 ³ copies/ml
F1-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Candida albicans Real-TM Real Time Amplification kit	R	€€	100	1 x10 ³ copies/ml
TF3-100FRT SA, RG, IQ, MX, SC, A, B	Candida albicans/C.glabrata/C.krusei Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	1 x10 ³ copies/ml
F3-100FRT SA, RG, IQ, MX, SC, A, B	Candida albicans/C.glabrata/C.krusei Real-TM Real Time Amplification kit	R	€€	100	1 x10 ³ copies/ml
TF5-100FRT SA, RG, B	Candidosis Real-TM Quant Complete Real Time Amplification kit for detection of <i>C.albicans</i> , <i>C.glabrata</i> , <i>C.krusei</i> , <i>C.parapsilosis</i> , <i>C. tropicalis</i> with DNA-Sorb-A extraction kit	R		100	1 x10 ³ copies/ml
F5-100FRT SA, RG, B	Candidosis Real-TM Quant Real Time Amplification kit for detection of <i>C.albicans</i> , <i>C.glabrata</i> , <i>C.krusei</i> , <i>C.parapsilosis</i> , <i>C. tropicalis</i>	R		100	1 x10 ³ copies/ml
F-1-100R	<i>Candida albicans</i> 500/730 IC	A		110	2,5 x10 ³ copies/ml

Treponema pallidum

Treponema pallidum is a species of spirochaete bacterium with subspecies that cause a disease such as syphilis. Syphilis is a sexually transmitted disease whose route of transmission is almost always through sexual contact, although there are examples of congenital syphilis via transmission from mother to child in utero or at birth. *T. pallidum* is transmitted via penetration of the spirochetes through mucosal membranes and abrasions on epithelial surfaces. Incubation time from exposure to development of primary lesions, which occur at the primary site of inoculation, averages 3 weeks but can range from 10-90 days. Syphilis can generally be treated with antibiotics, including penicillin. If left untreated, syphilis can damage the heart, aorta, brain, eyes, and bones. In some cases these effects can be fatal.



Treponema pallidum Kits

TB20-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Treponema pallidum Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	100	5 x10 ² copies/ml
B20-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Treponema pallidum Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
B-20-100R	<i>Treponema pallidum</i> 273/668 IC	A	€€	110	1 x10 ³ copies/ml

Bacterial vaginosis

Bacterial vaginosis (BV) is considered to be the most common cause of vaginal inflammation among both pregnant and non-pregnant women and prevalences between 4.9% and 36.0% have been reported from European and American studies. It previously was called nonspecific vaginitis or Gardnerella-associated vaginitis.

The adult human vagina is a complex ecosystem containing an abundance of microorganisms. In women of childbearing age this system is dominated by *Lactobacillus spp.*, a genus of gram-positive, nonmotile rod-like bacteria, a defining characteristic of which is the ability to grow in acid media and tolerate acid conditions (pH < 4.5); lactobacilli also ferment carbohydrates to produce lactic acid and produce H₂O₂ which provides a natural defense against *Gardnerella vaginalis*. In bacterial vaginosis (BV) the balance of flora is changed with reduced numbers of lactobacilli (normal concentration 10⁶ – 10¹⁰ CFU/ml) and an increase in numbers of other facultative and anaerobic species such as anaerobic cocci *Prevotella spp.*, *Gardnerella vaginalis*, and *Mobiluncus spp.* (normal concentration < 10³-10⁵ CFU/ml). *G. vaginalis* is virtually always present at high concentrations in women who have BV but is also detected frequently in normal women and in some cases the concentration of *Gardnerella vaginalis* can reach 10⁷-10⁸ CFU/ml also in absence of BV, so the most important maker of BV is the ratio of logarithm concentration *Lactobacillus spp* and *G. vaginalis*.

The clinical significance of studying vaginal flora is that it helps determine the quantity of microorganisms and assess the ratio between the different groups of conditionally pathogenic microorganisms and the normal flora. The total quantity of bacteria serves as an indicator of infection level in the vaginal environment: under normal conditions it can vary between 10⁶ and 10⁹ (6-9 Log). The ratio between lactobacilli and the total bacterial quantity can be used as an indicator of the balance between the normal and conditionally pathogenic flora: the normal proportion of lactobacilli should be 95 to 100% of the total bacterial quantity.

Bacterial vaginosis Kits

TR-B7-100FRT SA, RG, IQ, SC, MX, A, B	G.vaginalis/Lactobac. species Real-TM Quant Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	2,5 x10 ³ copies/ml
R-B7-100FRT SA, RG, IQ, SC, MX, A, B	G.vaginalis/Lactobac. species Real-TM Quant Real Time Amplification kit	R	€€	100	2,5 x10 ³ copies/ml
TB7-100FRT SA, RG, IQ, MX, SC, A, B	Gardnerella vaginalis Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	2,5 x10 ³ copies/ml
B7-100FRT SA, RG, IQ, MX, SC, A, B	Gardnerella vaginalis Real-TM Real Time Amplification kit	R	€€	100	2,5 x10 ³ copies/ml
TB74-100FRT SA, RG, IQ, MX, SC, A, B	Bacterial Vaginosis Real-TM Quant Complete Real Time Test for quantitative detection of Gardnerella vaginalis, Atopobium vaginae, Lactobacillus spp. and total bacteriae quantity in the vaginal biotope with DNA-Sorb-A extraction kit	R	€€	100	2,5 x10 ³ copies/ml
B74-100FRT SA, RG, IQ, MX, SC, A, B	Bacterial Vaginosis Real-TM Quant Multiplex RT-PCR for quantitative detection of Gardnerella vaginalis, Atopobium vaginae, Lactobacillus spp. and total bacteriae quantity in the vaginal biotope	R	€€	100	2,5 x10 ³ copies/ml
B-7-50R	Gardnerella vaginalis 355/640 IC	A		55	2,5 x10 ³ copies/ml
B-7-100R	Gardnerella vaginalis 355/640 IC	A	€€	110	2,5 x10 ³ copies/ml

Men's Health

Chronic bacterial prostatitis (CBP) represents a bacterial infection of the prostate gland. CBP causes an associated symptom complex, the hallmark of which is the occurrence of relapsing urinary tract infections. Approximately half of all men eventually develop symptoms consistent with prostatitis. This symptom complex accounts for approximately 25% of urologic evaluations in men.

Men's Health Kits

01765-50 SA, RG, IQ, MX, A, B	ProstateScreen Real-TM NEW Real Time PCR test for diagnosis of chronic bacterial prostatitis (<i>Enterobacter spp.</i> , <i>Klebsiella spp.</i> , <i>Enterococcus faecalis</i> , <i>E. faecium</i> , <i>Escherichia coli</i> , <i>Proteus spp.</i> , <i>Pseudomonas aeruginosa</i> , <i>Serratia spp.</i> , <i>Staphylococcus aureus</i> , <i>Streptococcus spp.</i>) 12 x 8 strip format	R		96	
----------------------------------	--	---	--	----	--

Cytomegalovirus (CMV)

Cytomegalovirus (CMV) is a double-stranded DNA virus and is a member of the Herpesviridae family. The other family members include herpes simplex virus type 1 (HSV-1 or HHV-1) and herpes simplex virus type 2 (HSV-2 or HHV-2), varicella zoster virus (VZV), human herpesvirus (HHV)–6, HHV-7, and HHV-8. CMV shares many attributes with other herpes viruses, including genome, virion structure, and the ability to cause latent and persistent infections. Human CMV grows only in human cells and replicates best in human fibroblasts. About 58.9% of individuals aged 6 and over are infected with CMV while 90.8% of individuals aged 80 and over are positive for CMV antibodies. Symptomatic CMV disease in immunocompromised individuals can affect almost every organ of the body, resulting in fever of unknown origin, pneumonia, hepatitis, encephalitis, myelitis, colitis, uveitis, retinitis, and neuropathy. In patients coinfecting with HIV, CMV infection leads to progression to AIDS and eventually death, even in those receiving highly active antiretroviral therapy (HAART).

CMV Kits

TV7-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	CMV Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R		100	5 x10 ² copies/ml
V7-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	CMV Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
TV7-100/2FRT SA, RG, IQ, SC, MX, A, IL, B	CMV Real-TM Quant Complete Real Time Test with DNA-Sorb-B extraction kit	R		100	2 x10 ² copies/ml
TV7-100/2FRT C SA, RG, IQ, SC, MX, A, IL, B	CMV Real-TM Quant Complete Real Time Test with with Column extraction kit	R		100	2 x10 ² copies/ml
V7-100/2FRT SA, RG, IQ, SC, MX, A, IL, B	CMV Real-TM Quant Real Time PCR kit with the DNA extraction controls	R	€€	100	2 x10 ² copies/ml
TV48-100FRT SA, RG, IQ, SC, MX, A	CMV/EBV/HHV6 Real-TM Quant Complete Real Time Test with DNA-Sorb-B extraction kit	R		100	2 x10 ² copies/ml
V48-100FRT SA, RG, IQ, SC, MX, A	CMV/EBV/HHV6 Real-TM Quant Real Time Amplification kit	R	€€	100	2 x10 ² copies/ml
V-7-100R	CMV 500/800 IC	A		110	5 x10 ² copies/ml

Herpes simplex

Herpes simplex is a viral disease caused by both herpes simplex virus type 1 (HSV-1) and type 2 (HSV-2). Infection with the herpes virus is categorized into one of several distinct disorders based on the site of infection. Oral herpes, the visible symptoms of which are colloquially called *cold sores* or *fever blisters*, infects the face and mouth. Oral herpes is the most common form of infection. Genital herpes, known simply as *herpes*, is the second most common form of herpes. Other disorders such as herpetic whitlow, herpes gladiatorum, ocular herpes (keratitis), cerebral herpes infection encephalitis, Mollaret's meningitis, neonatal herpes, and possibly Bell's palsy are all caused by herpes simplex viruses. Varicella zoster virus (VZV) is an herpes viruses known to infect humans (and other vertebrates) that commonly causes chicken-pox in children.

Herpes simplex Kits

TV8-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	HSV 1/2 Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	5 x10 ² copies/ml
V8-100FRT SA, RG, IQ, SC, MX, A, IL, B, LC	HSV 1/2 Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
TV38-100FRT SA, RG, IQ, MX, SC, A, IL, B	HSV 1/2 Typing Real-TM Complete Real Time Test with DNA-Sorb-A extraction kit	R	€€	100	5 x10 ² copies/ml

Herpes simplex Kits

V38-100FRT SA, RG, IQ, SC,MX, A,IL,B	HSV 1/2 Typing Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
TV60-100FRT SA, RG, IQ, SC,MX, A,IL,B	HSV/CMV Screen Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R		100	5 x10 ² copies/ml
V60-100FRT SA, RG, IQ, SC,MX, A,IL,B	HSV/CMV Screen Real-TM Real Time Amplification kit	R	€€	100	5 x10 ² copies/ml
TV61-50FRT SA, RG, IQ, SC,MX, A,IL,B	VZV Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	5 x10 ² copies/ml
V61-50FRT SA, RG, IQ, SC,MX, A,IL,B	VZV Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
V-8-100R	HSV 430/720 IC	A	€€	110	1 x10 ³ copies/ml

Epstein-Barr virus (EBV)

Epstein-Barr virus (EBV), or human herpesvirus 4, is a gammaherpesvirus that infects more than 95% of the world's population. The most common manifestation of primary infection with this organism is acute infectious mononucleosis, a self-limited clinical syndrome that most frequently affects adolescents and young adults. Classic symptoms include sore throat, fever, and lymphadenopathy. However, Epstein-Barr virus is also a human tumor virus, the first virus associated with human malignancy (nasopharyngeal carcinoma and Burkitt lymphoma).

EBV Kits

TV9-50FRT SA, RG, IQ, SC,MX, A,IL,B	EBV Real-TM Quant Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	2 x10 ² copies/ml
V9-50FRT SA, RG, IQ, SC,MX, A,IL,B	EBV Real-TM Quant Real Time Amplification kit with the DNA extraction controls	R	€€	50	2 x10 ² copies/ml
V-9-100R	EBV 290	A	€€	110	5 x10 ² copies/ml

Human Herpes Virus 6 (HHV6)

Human herpesvirus 6 (HHV-6) was the sixth herpesvirus discovered. Isolated in 1986 during attempts to find novel viruses in patients with lymphoproliferative diseases, HHV-6 is now recognized as a T-cell lymphotropic virus with high affinity for CD4 lymphocytes. A beta herpesvirus (like cytomegalovirus [CMV] and human herpesvirus type 7), HHV-6 has two variants, A and B. HHV-6B causes the childhood illness roseola infantum, while HHV-6A has been isolated mainly in immunocompromised hosts. HHV-6 also has been implicated in the pathogenesis of white-matter demyelination in persons with AIDS dementia complex; however, causality has not been proven. HHV-6 has been isolated from various tissues, cells, and fluid in association with the following conditions: Kikuchi lymphadenitis, Lymphoma, Lymphadenopathy, Sjögren syndrome, Sarcoidosis, Systemic lupus erythematosus, Guillain-Barré syndrome, Multiple sclerosis.

HHV-6

TV10-100FRT SA, RG, IQ, SC,MX, A,IL,B	HHV6 Real-TM Quant Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	100	2 x10 ² copies/ml
V10-100FRT SA, RG, IQ, SC,MX, A,IL,B	HHV6 Real-TM Quant Real Time Amplification kit	R	€€	100	2 x10 ² copies/ml
V-10-100R	HHV 6 - 380	A	€€	110	5 x10 ² copies/ml

Rubella

The name rubella is derived from a Latin term meaning “little red.” Rubella is generally a benign communicable exanthematous disease. It is caused by **Rubella virus**, which is a member of the Rubivirus genus of the family Togaviridae. Nearly one half of individuals infected with this virus are asymptomatic. Infection in younger children is characterized by mild constitutional symptoms, rash, and suboccipital adenopathy; conversely, in older children, adolescents, and adults, rubella may be complicated by arthralgia, arthritis, and thrombocytopenic purpura. Rare cases of rubella encephalitis have also been described in children. The major complication of rubella is its teratogenic effects when pregnant women contract the disease, especially in the early weeks of gestation. The virus can be transmitted to the fetus through the placenta and is capable of causing serious congenital defects, abortions, and stillbirths.

Rubella virus Kits

TV24-50FRT SA, RG, IQ, SC, MX, A, IL, B	Rubella Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R		50	1 x10 ³ copies/ml
V24-50FRT SA, RG, IQ, SC, MX, A, IL, B	Rubella Real-TM Real Time Amplification kit	R	€	50	1 x10 ³ copies/ml
V-24-50R	Rosolia 302 (RT-PCR)	A		55	

Toxoplasma gondii

Toxoplasmosis is caused by infection with **Toxoplasma gondii**, an obligate intracellular parasite. The infection produces a wide range of clinical syndromes in humans, land and sea mammals, and various bird species. In most immunocompetent individuals, primary or chronic (latent) *T gondii* infection is asymptomatic. A small percentage of these patients eventually develop chorioretinitis, lymphadenitis, or, rarely, myocarditis and polymyositis. However, certain individuals are at high risk for severe or life-threatening toxoplasmosis. Individuals at risk for toxoplasmosis include fetuses, newborns, and immunologically impaired patients. Congenital toxoplasmosis is usually a subclinical infection. Among immunodeficient individuals, toxoplasmosis most often occurs in those with defects of T-cell-mediated immunity, such as those with hematologic malignancies, bone marrow and solid organ transplants, or AIDS.

Toxoplasma gondii Kits

TP1-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Toxoplasma gondii Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R		50	4 x10 ² copies/ml
P1-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Toxoplasma gondii Real-TM Real Time Amplification kit	R	€	50	4 x10 ² copies/ml

Parvovirus B19

Parvovirus B19 (B19V) is a single-stranded DNA virus of the family Parvoviridae and genus Erythrovirus. Human parvovirus B19 was shown to be the etiologic agent of erythema infectiosum in hematologically normal persons. Erythema infectiosum was originally named Fifth disease because it was the fifth of 6 classic exanthematous diseases of childhood to be described. Later, cases of nonimmune hydrops fetalis were reported when infection in a woman occurred during pregnancy. Patients who are immunocompromised (eg, receiving chemotherapy or immunosuppressive drugs or have immune defects [congenital and acquired]) may develop chronic parvovirus B19 infection that results in chronic anemia. Pure red cell aplasia (PRAC) persists until the virus is cleared and should be distinguished from the transient anemia described above. Chronic parvovirus B19 infection in transplant recipients has been linked to anemia, other hematologic abnormalities, myocarditis, and pneumonitis.

Parvovirus B19 Kits

TV49-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Parvovirus B19 Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R	€	50	2 x10 ² copies/ml
V49-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Parvovirus B19 Real-TM Real Time Amplification kit	R	€	50	2 x10 ² copies/ml

Influenza Virus

Influenza virus infection, one of the most common infectious diseases, is a highly contagious airborne disease that causes an acute febrile illness and results in variable degrees of systemic symptoms, ranging from mild fatigue to respiratory failure and death. These symptoms contribute to significant loss of workdays, human suffering, mortality, and significant morbidity. Influenza results from infection with 1 of 3 basic types of influenza virus—A, B, or C—which are classified within the family Orthomyxoviridae. These single-stranded RNA viruses are structurally and biologically similar but vary antigenically. The most common prevailing influenza A subtypes that infect humans are H1N1 and H3N2.



In 1997, an Avian subtype of influenza A, H5N1, was first described in Hong Kong. In 2008, more than 390 human cases had been documented and more than 246 persons had died following H5N1 outbreaks among poultry and resulting bird-to-human transmission.

On April 26, 2009, the US Department of Health and Human Services issued a nationwide public health emergency regarding Swine Influenza A (H1N1) virus infections in humans. As of early June 2009, H1N1 influenza had infected 28,774 people in 74 countries, and 144 deaths were confirmed to have been caused by the disease.

Influenza Kits

TV36-50FRT SA, RG, IQ, SC, MX, A, IL, B	Influenza A, B Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R	CE	50	1 x10 ³ copies/ml
V36-50FRT SA, RG, IQ, SC, MX, A, IL, B	Influenza A, B Real-TM Real Time Amplification kit	R	CE	50	1 x10 ³ copies/ml
TR-V33-FRT SA, RG, IQ, SC, MX, A, IL, B	Avian A Screening & Avian H5N1 Typing FRT Complete Real Time Test with Ribo-Sorb extraction kit	R	CE	50	1 x10 ³ copies/ml
R-V33-FRT SA, RG, IQ, SC, MX, A, IL, B	Avian A Screening & Avian H5N1 Typing FRT Real Time Amplification kit	R	CE	50	1 x10 ³ copies/ml
TV31-50FRT SA, RG, IQ, B	Avian A Screening & Avian H5 H7 Typing FRT Complete Real Time Test with Ribo-Sorb extraction kit	R		50	1 x10 ³ copies/ml
V31-50FRT SA, RG, IQ, B	Avian A Screening & Avian H5 H7 Typing FRT Real Time Amplification kit	R		50	1 x10 ³ copies/ml
TV55-50FRT SA, RG, IQ, SC, MX, A, IL, B	Swine Influenza Virus H1 Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R	CE	50	5 x10 ² copies/ml
V55-50FRT SA, RG, IQ, SC, MX, A, IL, B	Swine Influenza Virus H1 Real-TM Real Time Amplification kit	R	CE	50	5 x10 ² copies/ml
TV47-50FRT SA, RG, IQ, SC, MX, A, IL, B	Influenza A H5 H7 H9 Typing FRT Complete Real Time Test with Ribo-Sorb extraction kit	R	CE	50	1 x10 ³ copies/ml
V47-50FRT SA, RG, IQ, SC, MX, A, IL, B	Influenza A H5 H7 H9 Typing FRT Real Time Amplification kit	R	CE	50	1 x10 ³ copies/ml
TV54-50FRT SA, RG, IQ, B	Influenza A H1N1 & H3N2 Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R		50	1 x10 ³ copies/ml
V54-50FRT SA, RG, IQ, B	Influenza A H1N1 & H3N2 Real-TM Real Time Amplification kit	R		50	1 x10 ³ copies/ml
TV51-50FRT RG, IQ, SA, B	Parainfluenza Virus Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R	CE	50	1 x10 ³ copies/ml
V51-50FRT RG, IQ, SA, B	Parainfluenza Virus Real-TM Real Time Amplification kit	R	CE	50	1 x10 ³ copies/ml

Chlamydia & Mycoplasma pneumoniae

Chlamydophila (formerly *Chlamydia*) *pneumoniae* causes mild pneumonia or bronchitis in adolescents and young adults. Older adults may experience more severe disease and repeated infections. Approximately 50% of young adults and 75% of elderly persons have serological evidence of previous infection. The pathogen is estimated to cause 10-20% of community-acquired pneumonia cases among adults. The estimated number of cases of C pneumoniae pneumonia is 300,000 cases per year.

Chlamydia pneumoniae Kits

TB42-4-50FRT SA, RG, IQ, B	Mycoplasma pneumoniae / Chl. pneumoniae Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R		50	5 x10 ² copies/ml
B42-4-50FRT SA, RG, IQ, B	Mycoplasma pneumoniae / Chl. pneumoniae Real-TM Real Time Amplification kit	R		50	5 x10 ² copies/ml

Pseudomonas

Pseudomonas is a gram-negative rod that belongs to the family Pseudomonadaceae. More than half of all clinical isolates produce the blue-green pigment pyocyanin. These pathogens are widespread in nature, inhabiting soil, water, plants, and animals (including humans). *Pseudomonas aeruginosa* has become an important cause of infection, especially in patients with compromised host defense mechanisms. Pseudomonal bacteremia occurs in association with malignancy, chemotherapy, AIDS, burn wound sepsis, and diabetes. It is a frequent cause of nosocomial infections such as pneumonia, endocarditis, meningitis, urinary tract infections (UTIs) and bacteremia.

Pseudomonas aeruginosa Kits

B76-50FRT SA, RG, IQ, SC, MX, A, IL, B	Pseudomonas aeruginosa Real-TM Quant Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
---	--	---	----	----	------------------------------

Legionella pneumophila

Legionella pneumophila is a thin, pleomorphic, flagellated Gram-negative bacterium of the genus Legionella. L. pneumophila is the primary human pathogenic bacterium in this group and is the causative agent of legionellosis or Legionnaires' disease. Legionella pneumophila (named in memory of the deceased veterans) is ubiquitous to aquatic environments worldwide and resided as an intracellular parasite of amoeba and protozoa provided a link between natural environment and human disease. Thus, environmental monitoring, especially of potable water, cooling towers, and related sources, is a major focus in efforts to control the spread of this disease.

Since the initial identification of 235 cases in 1976, Legionnaires disease has become recognized as the most common cause of atypical pneumonia in hospitalized patients. It is the second most common cause of community-acquired bacterial pneumonia with 25% mortality rate.

Legionella pneumophila Kits

TB50-50FRT RG, SA, IQ, B	Legionella pneumophila Real-TM Real Time Amplification kit with DNA-Sorb-B extraction kit	R	€€	50	5 x10 ² copies/ml
B50-50FRT RG, SA, IQ, B	Legionella pneumophila Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml




Figure. Incidence of Legionella related to sex and age.

Human respiratory syncytial virus

Human respiratory syncytial virus (hRSV) is a negative-sense, single-stranded RNA virus of the family *Paramyxoviridae*, the most common cause of severe lower respiratory tract disease among infants and young children, typically infects persons by age 2 years and can cause subsequent infections throughout life. RSV infection primarily manifests as bronchiolitis or pneumonia and results in approximately 75,000 to 125,000 hospitalizations in the United States each year. Persons at increased risk for severe disease or death include premature infants, older adults, and persons of any age with compromised respiratory, cardiac, or immune systems.

hRSV Kits


TV37-50FRT RG	hRSV Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R	€€	50	5 x10 ² copies/ml
V37-50FRT RG	hRSV Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
V-37-50R	Respiratory syncytial virus (hRSV) 298/550 IC (RT-PCR)	A	€€	55	5 x10 ² copies/ml

Acute Respiratory Viral Infections

ARVI Screen Real-TM

Real Time PCR kit for detection of 11 respiratory pathogens in one clinical sample:

- human parainfluenza virus-1-4 (hPiv) RNA;
- HKUI human coronavirus (hCov) RNA;
- human rhinovirus (hRv) RNA;
- human B, C, and E adenovirus (hAdv) DNA;
- human bocavirus (hBov) DNA



ARVI Kits

TV57-100FRT SA, RG, IQ, MX, SC, A, B	ARVI Screen Real-TM Multiplex RT-PCR detection and identification of human parainfluenza virus-1-4 (hPiv) RNA; OC43, E229, NL63, and HKUI human coronavirus (hCov) RNA; human rhinovirus (hRv) RNA; human B, C, and E adenovirus (hAdv) DNA; and human bocavirus (hBov) DNA with Ribo-Sorb extraction kit	R	€€	100	
V57-100FRT SA, RG, IQ, MX, SC, A, B	ARVI Screen Real-TM Real Time Amplification kit	R	€€	100	

Coronavirus (MERS+SARS)

Coronaviruses are species in the genera of virus belonging to the subfamily Coronavirinae in the family Coronaviridae. Coronaviruses are enveloped viruses with a positive-sense RNA genome and with a nucleocapsid of helical symmetry. The most known human coronavirus, SARS-CoV which causes SARS, has a unique pathogenesis because it causes both upper and lower respiratory tract infections and can also cause gastroenteritis.

Coronavirus Kits

V65-50FRT SA, RG, IQ, MX, SC, A, B	MERS-CoV Real-TM Real Time Amplification kit for detection and differentiation of MERS-CoV (Middle East Respiratory Syndrome Coronavirus) and SARS-CoV (Severe Acute Respiratory Syndrome Coronavirus)	R	€€	50	
---------------------------------------	--	---	----	----	--

Bordetella pertussis

Pertussis, also named whooping cough, is a highly contagious bacterial disease caused by *Bordetella pertussis*. Symptoms are initially mild, and then develop into severe coughing, which produce the namesake high-pitched “whoop” sound in infected babies and children when they inhale air after coughing. The coughing stage lasts for approximately six weeks before subsiding.

Pertussis Kits


B84-100FRT SA, RG, MX, SC, IQ, A, B	Bordetella pertussis/B.parapertussis/B.bronchiseptica Real-TM Real Time Amplification kit	R	CE	100	1 x10 ³ copies/ml
--	---	---	----	-----	------------------------------

Tuberculosis

Tuberculosis (abbreviated as TB for tubercle bacillus) is a common and deadly infectious disease caused by mycobacteria, mainly *Mycobacterium tuberculosis*. Tuberculosis most commonly attacks the lungs (as pulmonary TB) but can also affect the central nervous system, the lymphatic system, the circulatory system, the genitourinary system, bones, joints and even the skin. Other mycobacteria such as *Mycobacterium bovis*, *Mycobacterium africanum* and *Mycobacterium microti* can also cause tuberculosis.

Early diagnosis of tuberculosis makes effective treatment possible and increases the probability of clinical outcome owing to quite effective antituberculosis therapy, however the tuberculosis diagnosis has certain difficulties. The application of molecular biology methods allow to overcome the difficulties in the diagnosis of *Mycobacterium tuberculosis*, but due to the biological peculiarities of this microorganism and immune response of human organism, tuberculosis cannot be diagnosed only by one method.

The development of test to differentiate between infection with *Mycobacterium tuberculosis* or *Mycobacterium bovis* and vaccination with *M. bovis* BCG could greatly assist in the diagnosis of early infection as well as enhance the use of tuberculosis vaccines on a wider scale.




TB Kits

TB15-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	MTB Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	CE	50	2 x10 ² copies/ml
B15-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	MTB Real-TM Real Time Amplification kit	R	CE	50	2 x10 ² copies/ml
TB41-50FRT SA, RG, IQ, SC, MX, A, IL, B	Mycobacterium tuberculosis Diff Real-TM Complete Test with Magno-Sorb-Tub extraction kit	R	CE	50	5 x10 ² copies/ml
B41-50FRT SA, RG, IQ, SC, MX, A, IL, B	Mycobacterium tuberculosis Diff Real-TM Real Time Amplification kit	R	CE	50	5 x10 ² copies/ml
B-15-50R	MTB complex 390/750 IC	A		55	5 x10 ² copies/ml
B-15-100R	MTB complex 390/750 IC	A	CE	110	5 x10 ² copies/ml
B-38-50R	Mycobacterium paratuberculosis 209/700 IC	A		55	5 x10 ² copies/ml

Enteroviruses

Human **Enteroviruses** are ubiquitous viruses transmitted from person to person via direct contact with virus shed from the gastrointestinal or upper respiratory tract. Enteroviruses belong to the Picornaviridae family of viruses and are traditionally divided into 5 subgenera based on differences in host range and pathogenic potential. Each subgenus contains a number of unique serotypes, which are distinguished basing on neutralization by specific antisera. The subgenera include polioviruses, coxsackievirus (groups A and B), and echoviruses.

Enteroviruses cause a wide range of infections. **Poliovirus**, the prototypical enterovirus, can cause a subclinical or mild illness, aseptic meningitis, or paralytic poliomyelitis, a disease that has been eradicated in the United States and other developed countries. The nonpolio viruses (group A and B coxsackieviruses, echoviruses, enteroviruses) continue to be responsible for a wide spectrum of diseases in persons of all ages, although infection and illness occur most commonly in infants.




Enteroviruses Kits

TV16-50FRT RG, iQ, SA, B	Enterovirus Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R	€	50	1 x10 ³ copies/ml
V16-50FRT RG, iQ, SA, B	Enterovirus Real-TM Real Time Amplification kit	R	€	50	1 x10 ³ copies/ml
V64-50FRT SA, RG, iQ, SC, MX, A, IL, B	Enterovirus 71-Type Real-TM Real Time Amplification kit for detection of Enterovirus 71	R	€	50	
TV58-50FRT SA, RG, iQ, SC, MX, A, IL, B	Poliovirus Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R	€	50	1 x10 ³ copies/ml
V58-50FRT SA, RG, iQ, SC, MX, A, IL, B	Poliovirus Real-TM Real Time Amplification kit	R	€	50	1 x10 ³ copies/ml
V-16-50R	Enterovirus 207 (RT-PCR)	A	€	55	1 x10 ³ copies/ml

Meningitis


Meningitis is a clinical syndrome characterized by inflammation of the meninges. Clinically, this medical condition manifests with meningeal symptoms (eg, headache, nuchal rigidity, photophobia) and an increased number of white blood cells in the cerebrospinal fluid (CSF). Depending on the duration of symptoms, meningitis may be classified as acute or chronic. Acute bacterial meningitis denotes a bacterial cause of this syndrome. Depending on the specific bacterial cause, the syndrome may be called, for example, Streptococcus pneumoniae meningitis, Neisseria meningitidis, or Haemophilus influenzae meningitis. Kit NHS Meningitis Real-TM is a Real-Time test for the detection and differentiation of Neisseria meningitidis, Haemophilus influenzae and Streptococcus pneumoniae in the biological materials. DNA is extracted from specimens, amplified using RT-amplification and detected using fluorescent reporter dye probes specific for N.meningitidis, H.influenzae, S.pneumoniae DNA and IC (Internal Control).



N.meningitidis



H.influenzae



Str. pneumoniae

Meningitis Kits

TB25-50FRT SA, RG, IQ, SC, MX, A, B	NHS Meningitidis (N.meningitidis, H.influenzae, Str.pneumoniae) Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€	50	5 x10 ² copies/ml
B25-50FRT SA, RG, IQ, SC, MX, A, B	NHS Meningitidis (N.meningitidis, H.influenzae, Str.pneumoniae) Real-TM Real Time Amplification kit	R	€	50	5 x10 ² copies/ml
B-14-100FRT SA, RG, IQ, SC, MX, A, B	Listeria monocytogenes Real-TM Quant Real Time Amplification kit	R	€	100	1 x10 ³ copies/ml
B-25-50R	NHS Meningitidis (N.meningitidis, H.influenzae, Str.pneumoniae)	A	€	55	1 x10 ³ copies/ml
B-26-50R	Nesseria meningitidis A, B, C types	A		55	1 x10 ³ copies/ml
B-14-50R	Listeria monocytogenes 280	A		55	1 x10 ³ copies/ml

Intestinal Infections

Acute Intestinal Infections (A.I.I) are one of the primary causes of hospitalization in infectious disease departments. In accordance with the data provided by the contemporary literature the following bacterial and viral agents are the most often detectable and generally spread etiological agents of All:

1. Bacterial agents:
 - Shigella species microorganisms and enteroinvasive E coli (EIEC);
 - Salmonella species microorganisms;
 - Thermophilic group of Campylobacter species microorganisms;
 - Enteropathogenic E coli (EPEC) and enteroaggregative E coli (EAEC);
2. Viral agents
 - Group A rotaviruses;
 - Genotype 2 noroviruses;
 - Group F adenoviruses (Types 40 and 41);
 - Astroviruses.

The following causative agents are less widely or not universally spread but are no less important for epidemic outbreaks:

1. Vibrio cholerae;
2. Yersinia pseudotuberculosis;
3. Clostridium difficile;
4. Enterotoxigenic E. coli (ETEC), Enterohemorrhagic E. coli (EHEC);
5. Genotype 1 Enteroviruses;
6. Group C Rotaviruses.

Intestinal Infections Kits


TV40-50FRT SA, RG, IQ, SC, MX, A, B	Rotavirus/Norovirus/Astrovirus Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R	€€	50	5 x10 ² copies/ml
V40-50FRT SA, RG, IQ, SC, MX, A, B	Rotavirus/Norovirus/Astrovirus Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
TB44-50FRT SA, RG, IQ, SC, MX, A, B	Shigella/Salmonella/Campylobacter Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	5 x10 ² copies/ml
B44-50FRT SA, RG, IQ, SC, MX, A, B	Shigella/Salmonella/Campylobacter Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
TB45-50FRT SA, RG, IQ, SC, MX, A, B	A.I.I. (Acute Intestinal Infections) Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R	€€	50	5 x10 ² copies/ml
B45-50FRT SA, RG, IQ, SC, MX, A, B	A.I.I. (Acute Intestinal Infections) Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
TB62-50FRT SA, RG, IQ, MX, SC, A, B	Escherichioses Screen & Diff Real-TM Multiplex RT-PCR detection and identification E.coli (EPEC, ETEC, EIEC, EHEC, and EA _g EC)	R	€€	50	5 x10 ² copies/ml
B62-50FRT SA, RG, IQ, MX, SC, A, B	Escherichioses Screen & Diff Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
B64-50FRT SA, RG, IQ, SC, MX, A, B	Yersinia enterocolitica/Y.pseudotuberculosis Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
B88-100FRT SA, RG, IQ, SC, MX, A, B	Aerobic complex Real-TM Real Time Amplification kit for detection of Enterobacteriaceae (E. coli, Klebsiella spp, Proteus spp, Streptococcus spp)	R	€€	100	
V-15-50R	Rotavirus 280 (RT-PCR)	A	€€	55	1 x10 ³ copies/ml
V-19-50R	Astrovirus 175 (RT-PCR)	A	€€	55	1 x10 ³ copies/ml
B-23-50R	Clostridium difficile 420	A		55	1 x10 ³ copies/ml
B-39-50R	Yersinia pseudotuberculosis 527	A		55	1 x10 ³ copies/ml
B-22-50R	Yersinia Enterocolitica 270	A		55	1 x10 ³ copies/ml

Intestinal Infections Kits: Food Pathogens

TB35-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Campylobacter species Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	5 x10 ² copies/ml
B35-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Campylobacter species Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
TB58-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Cronobacter sakazakii Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	5 x10 ² copies/ml
B58-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Cronobacter sakazakii Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
TB12-50FRT RG, IQ, SC, MX, A, B	Shigella Spp & EIEC Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	5 x10 ² copies/ml
B12-50FRT SA, RG, IQ, SC, MX, A, B	Shigella Spp & EIEC Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
TB59-50FRT SA, RG, IQ, SC, MX, A, B	EHEC Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	5 x10 ² copies/ml
B59-50FRT SA, RG, IQ, SC, MX, A, B	EHEC Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
TB11-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Salmonella species Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	5 x10 ² copies/ml
B11-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	Salmonella species Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml
B63-50FRT SA, RG, B, SC	Salmonella Diff Real-TM Real Time Amplification kit for detection and differentiation of <i>S. typhi</i> , <i>S. paratyphi</i> , <i>Salmonella</i> H1 flagellar antigen types (<i>S. stanley</i> , <i>S. isangi</i> , <i>S. muenchen</i> , <i>S. gaminara</i> , <i>S. utrecht</i>)	R		50	
B-11-50R	Salmonella species 225/670 IC	A		55	1 x10 ³ copies/ml
B-12-50R	Shigella species 400/800 IC	A		55	1 x10 ³ copies/ml
B-13-50R	Campylobacter jejuni 415	A		55	1 x10 ³ copies/ml
B-35-50R	Campylobacter species 520	A		55	1 x10 ³ copies/ml

Helicobacter Pylori

Helicobacter pylori is a gram-negative bacillus responsible for one of the most common infections found in humans worldwide. H pylori organisms are spiral-shaped gram-negative bacteria that are highly motile because of multiple unipolar flagella. They are microaerophilic and potent producers of the enzyme urease. H pylori inhabits the mucus adjacent to the gastric mucosa.




Helicobacter Pylori Kits

TB9-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	H.pylori Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€	50	1 x10 ³ copies/ml
B9-50FRT SA, RG, IQ, SC, MX, A, IL, B, LC	H.pylori Real-TM Real Time Amplification kit	R	€	50	1 x10 ³ copies/ml
B-9-50R	Helicobacter pylori 520	A	€	55	1 x10 ³ copies/ml

Tropical parasites

Parasites are one of the four main agents of infection known to man, the other three being bacterial, fungal and viral. Tropical parasites are prevalent in or unique to tropical and subtropical regions. The diseases are less prevalent in temperate climates, partly because of the occurrence of a cold season, which controls the insect population by forcing hibernation. Insects such as mosquitoes and flies are one of the most common disease carrier. These insects may carry a parasite that is infectious to humans and animals.



Tropical Parasites Kits

N3-50FRT SA, RG, IQ, SC, MX, A	Leishmania spp. Real-TM Real Time Amplification kit	R		50	1 x10 ³ copies/ml
-----------------------------------	---	---	--	----	------------------------------

Lyme disease


Lyme disease (LD) is a vector-borne, multisystem inflammatory disease caused by the spirochete ***Borrelia burgdorferi sensu lato***. It is transmitted to humans by infected ticks of the *Ixodes* genus. After entering the circulation, the organism invades the cutaneous, synovial, cardiac, and nervous systems. Spirochetes have also been demonstrated histologically in bone marrow, the spleen, lymph nodes, the liver, testes, and the placenta during early hematogenous dissemination.

Borrelia burgdorferi Kits

TB37-50FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Borrelia burgdorferi sensu lato Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	1 x10 ³ copies/ml
B37-50FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Borrelia burgdorferi sensu lato Real-TM Real Time Amplification kit	R	€€	50	1 x10 ³ copies/ml
B-37-50R	Borrelia burgdorferi sensu lato	A	€€	55	

Brucellosis

Brucellosis is a worldwide zoonosis caused by infection with the bacterial genus *Brucella*. These organisms, which are small aerobic intracellular coccobacilli, localize in the reproductive organs of host animals, causing abortions and sterility. They are shed in large numbers in the animal's urine, milk, placental fluid, and other fluids. Exposure to infected animals and animal products causes brucellosis in humans. The global burden of human brucellosis remains enormous; it causes more than 500,000 infections per year worldwide. Among the 4 *Brucella* species known to cause disease in humans (*B abortus*, *B melitensis*, *B canis*, *B suis*), *B melitensis* is thought to be the most virulent and causes the most severe and acute cases of brucellosis.



Brucella Kits

TB10-50FRT SA, RG, IQ, MX, B	Brucella Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	1 x10 ³ copies/ml
B10-50FRT SA, RG, IQ, MX, B	Brucella Real-TM Real Time Amplification kit	R	€€	50	1 x10 ³ copies/ml
B-10-50R	Brucella species 460/770 IC	A	€€	55	1 x10 ³ copies/ml

Dengue virus

Dengue virus (DENV) is the cause of dengue fever. It is a mosquito-borne single positive-stranded RNA virus of the family Flaviviridae; genus Flavivirus. Four serotypes can cause the full spectrum of disease. There are not yet any vaccines to prevent infection with dengue virus and the most effective protective measures are those that avoid mosquito bites. When infected, early recognition and prompt supportive treatment can substantially lower the risk of developing severe disease.

Dengue Kits

V63-G-50FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Dengue genotype Real-TM Real Time PCR kit for detection and differentiation of Dengue genotypes 1, 2, 3 and 4	R	€€	50	500 copies/ml
V63-S-50FRT	Dengue Real-TM Real Time PCR kit for detection of Dengue Virus	R	€€	50	500 copies/ml

Pestis

Yersinia pestis is a Gram-negative rod-shaped bacterium. It is a facultative anaerobe able to infect humans and other animals. Many evidence suggest that it was a contributing factor in many plagues throughout human history. The reservoir commonly associated with *Y. pestis* is several species of rodents (marmot, rats).

Pestis Kits

B79-50FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Yersinia pestis Real-TM Real Time Amplification kit	R	€€	50	5 x10 ³ copies/ml
--	---	---	----	----	------------------------------

Bacillus Anthracis

Bacillus anthracis is a Gram-positive spore-forming, rod-shaped bacterium, with a width of 1-1.2µm and a length of 3-5µm. Three forms of anthrax disease are recognized based on their form of inoculation: Cutaneous - the most common form (95%), causes a localized inflammatory black necrotic lesion (eschar), Pulmonary - highly fatal and characterized by sudden massive chest edema followed by cardiovascular shock, Gastrointestinal - rare but also fatal (causes death to 25%) type results from ingestion of spores

Bacillus Anthracis Kits

TB101-50FRT SA, RG, IQ, MX, B	Bacillus anthracis Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	5 x10 ² copies/ml
B101-50FRT SA, RG, IQ, MX, B	Bacillus anthracis Real-TM Real Time Amplification kit	R	€€	50	5 x10 ² copies/ml

Vibrio Cholerae

Cholera is caused by *Vibrio cholerae*, the most feared epidemic diarrheal disease because of its severity. Dehydration and death can occur within hours of infection. Robert Koch discovered *V cholerae* in 1883 during an outbreak in Egypt. The organism is a comma-shaped, gram-negative aerobic bacillus whose size varies from 1-3 µm in length by 0.5-0.8 µm in diameter. Its antigenic structure consists of a flagellar H antigen and a somatic O antigen.

Cholera Kits

TB53-50FRT SA, RG, IQ, MX, B	Vibrio cholerae Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R		50	5 x10 ² copies/ml
B53-50FRT SA, RG, IQ, MX, B	Vibrio cholerae Real-TM Real Time Amplification kit	R		50	5 x10 ² copies/ml

Leptospira

Leptospirosis is a worldwide zoonosis caused by pathogenic species of the genus *Leptospira*. In 90% of cases, leptospirosis manifests as an acute febrile illness with a biphasic course and an excellent prognosis. Nonspecific signs and symptoms of leptospirosis (eg, fever, headache, nausea, vomiting) are often confused with viral illness. In 10% of cases, the presentation is more dramatic, and the infection has a mortality rate of 10%. Known as Weil disease or icteric leptospirosis, the classic definition of this form of leptospirosis includes fever, jaundice, renal failure, and hemorrhage.

Leptospirosis Kits

TB49-50FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Leptospira 16s RNA Real-TM Complete Real Time Test with DNA-Sorb-B extraction kit	R	€€	50	1 x10 ³ copies/ml
B49-50FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Leptospira 16s RNA Real-TM Real Time Amplification kit	R	€€	50	1 x10 ³ copies/ml

Congo Crimea

Congo-Crimea hemorrhagic fever is a widespread tick-borne viral disease, a zoonosis of animals that may affect humans. The pathogenic virus, commonly present in East and West Africa, is a member of the Bunyaviridae family of RNA viruses. Clinical disease is rare in infected mammals, but commonly severe in infected humans, with a 30% mortality rate.

Congo Crimea Kits

V22-50FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Congo Crimea Real-TM Real Time Amplification kit	R	€€	50	1 x10 ³ copies/ml
--	--	---	----	----	------------------------------

Tick-Borne diseases

Ticks are arachnids, relatives of spiders that commonly live in wooded areas, brushy fields. They survive by eating blood from their hosts and they can pass infections from one host to another, including humans. Common tick-borne diseases are: **Tick-Borne Encephalitis (TBE), Anaplasmosis, Ehrlichiosis.**

TBEV Kits

V59-100FRT SA, RG, IQ, SC, MX, A, B	TBEV, B.burgdorferi, A.phagocytophilum, E.chaffeensis / E.muris Real-TM Real Time Amplification kit	R	€€	100	1 x10 ³ copies/ml
--	---	---	----	-----	------------------------------

West Nile Virus

The West Nile virus is one of the many members of the genus *Flavivirus* that are known to cause human disease. The life cycle of the West Nile virus involves the microbe's transmission from nonhuman animals to humans by way of *Aedes*, *Culex*, or *Anopheles* mosquitoes. The West Nile virus can infect horses, birds, dogs, and other mammals. The West Nile virus causes serious manifestations in approximately 1% of persons who are infected, with increased morbidity and mortality in individuals older than 50 years. In hospitalized patients, neurologic sequelae of the West Nile virus included severe muscle weakness, with approximately 10% of patients developing a complete flaccid paralysis. One in 150 West Nile virus infections results in encephalitis or meningitis, and the mortality rate from severe illness is 3-15%.

West Nile Virus Kits

TV53-50FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	West Nile Virus Real-TM Complete Real Time Test with Ribo-Sorb extraction kit	R	€€	50	5 x 10 ² copies/ml
V53-50FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	West Nile Virus Real-TM Real Time Amplification kit	R	€€	50	5 x 10 ² copies/ml

Coxiella burnetii

Q fever is a zoonosis caused by *Coxiella burnetii*, an obligate gram-negative intracellular bacterium. Most commonly reported in southern France and Australia, Q fever occurs worldwide. *C burnetii* infects various hosts, including humans, ruminants (cattle, sheep, goats), and pets. In rare cases, *C burnetii* infection in reptiles, birds, and ticks has been reported. *C burnetii* is excreted in urine, milk, feces, and birth products. These products, especially the latter, contain large numbers of bacteria that become aerosolized after drying. The bacterium is highly infectious, and only a few organisms can cause disease.

Coxiella burnetii Kits

B85-50FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Coxiella burnetii Real-TM Real Time Amplification kit	R	€€	50	1 x 10 ³ copies/ml
--	---	---	----	----	-------------------------------

Bacterial pyogenic infections and Pseudomonas aeruginosa

Pyogenic refers to bacterial infections that make pus, that is destroyed by bacteria such as *Streptococcus pyogenes* and *Staphylococcus aureus* through the release of leukocidins.

Bacterial pyogenic infections Kits

TB77-100FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Streptococcus B Real-TM Quant Complete Real Time Test with DNA/RNA Prep extraction kit	R	€€	100	
B77-100FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Streptococcus B Real-TM Quant Real Time Amplification kit	R	€€	100	
B82-100FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Streptococcus pyogenes Real-TM Quant Real Time Amplification kit	R	€€	100	
B76-50FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Pseudomonas aeruginosa Real-TM Quant Real Time Amplification kit	R	€€	50	

Ebola

Ebola virus disease (EVD), formerly known as Ebola haemorrhagic fever, is a severe, often fatal illness in humans. The virus is transmitted to people from wild animals and spreads in the human population through human-to-human transmission. The average EVD case fatality rate is around 50%. Case fatality rates have varied from 25% to 90% in past outbreaks.



Ebola Kits

V69-100FRT SA, RG, A, IQ, SC, MX, B, A, LC, IL	Ebola Zaire Real-TM Real Time PCR kit for detection of Ebola Zaire virus	R		100	
---	--	---	--	-----	--

Bacterial Drug Resistance

More commonly, the term Drug Resistance is used in the context of resistance that pathogens have “acquired”. When a microorganism is resistant to more than one drug, it is said to be multidrug-resistant (MDR). Resistance to β -lactam antibiotics which include the penicillins (oxacillin, methicillin, dicloxacillin, nafcillin etc.) and the cephalosporins make difficult to treat infections with standard types of antibiotics. In addition to β -lactam/carbapenem resistance, Enterobacteriaceae often carry genes that confer high levels of resistance to many other antimicrobials, often leaving very limited therapeutic options.

Drug resistance Kits

B78-100FRT SA, RG, MX, IQ, SC, A, B	MRSA Quant Real-TM Real Time Amplification kit	R	€€	100
C1-100FRT SA, RG, A, B	MDR MBL (VIM, IMP, NDM) Real-TM Real Time Amplification kit for detection of genes VIM, IMP and NDM in Enterobacteriaceae and NFGNB*	R	€€	100
C2-100FRT SA, RG, A, B	MDR KPC/OXA Real-TM Real Time Amplification kit for detection of genes KPC, OXA-48 and OXA-162 in Enterobacteriaceae and NFGNB	R	€€	100
T01781-50-T SA, RG, IQ, MX, A, B	Resistance cephalosporins-1 Real Time PCR test for Enterobacteriaceae resistance, genes CTX-M - ready to use 0,2 ml tube format	R		60
T01782-50-T SA, RG, IQ, MX, A, B	Resistance cephalosporins-2 Real Time PCR test for S.aureus resistance, genes Meca - ready to use 0,2 ml tube format	R		60
T01784-50-T SA, RG, IQ, MX, A, B	Resistance carbapenems-1 Real Time test for Enterobacteriaceae e Pseudomonas resistance, genes VIM - ready to use 0,2 ml tube format	R		60
T01785-50-T SA, RG, IQ, MX, A, B	Resistance carbapenems-2 Real Time PCR test for Enterobacteriaceae resistance, genes NDM - ready to use 0,2 ml tube format	R		60
T01786-50-T SA, RG, IQ, MX, A, B	Resistance carbapenems-3 Real Time PCR test for Enterobacteriaceae resistance, genes OXA-48 - ready to use 0,2 ml tube format	R		60
T01746-50-T SA, RG, IQ, MX, A, B	Resistance carbapenems-4 Real Time PCR test for Enterobacteriaceae resistance, genes KPC - ready to use 0,2 ml tube format	R		60
T01787-50-T SA, RG, IQ, MX, A, B	Resistance glycopeptides Real Time PCR test for Enterococcus faecalis and E. faecium resistance to Vancomycin and Teicoplanin, genes VanA and VanB - ready to use 0,2 ml tube format	R		60
T01788-50-T SA, RG, IQ, MX, A, B	Resistance penicillins Real Time PCR test for Enterobacteriaceae and Pseudomonas resistance to penicillins and first-generation cephalosporins, genes TEM - ready to use 0,2 ml tube format	R		60
T01747-50-T SA, RG, IQ, MX, A, B	Resistance macrolides-1 Real Time PCR test for Streptococcus spp Staphylococcus spp resistance to macrolides, lincosamides, streptogramins, gene ErmB - ready to use 0,2 ml tube format	R		60
T01748-50-T SA, RG, IQ, MX, A, B	Resistance macrolides-2 Real Time PCR test for Streptococcus spp resistance to macrolides; genes Mef - ready to use 0,2 ml tube format	R		60

*NFGNB: Non Fermenting Gram Negative Bacteria

Molecular Genetics

Cardiovascular diseases (CVD) are lifethreatening conditions which affect up to 10% of the human population. Thrombotic complications, such as an acute myocardial infarction, ischemic stroke, pulmonary embolism, deep venous thrombosis are the major causes of morbidity and mortality in the world. A wide spectrum of CVD with inherited genetic susceptibilities is now known and genetic susceptibility may be caused by mutations and single nucleotide polymorphisms in a variety of genes mainly involved in blood coagulation, regulation of blood pressure, and metabolism of lipids, glucose, homocysteine or iron.

Among the cardiovascular diseases markers have important role variations in the genes for blood coagulation factors V (**FV**), II (**protrombin**), XIII (**FXIII**), plasminogen activator inhibitor-1 (**PAI-1**), methylenetetrahydrofolate reductase (**MTHFR**), apolipoprotein B (**Apo B**), platelet glycoprotein IIIa (**GPIIIa**), β -fibrinogen (**FGB**) Moreover, an increased tendency to develop thrombosis, called also "thrombophilia", underlies the significant proportion of cases in the most common obstetric complications (recurrent pregnancy loss, fetal growth retardation, preeclampsia, abruptio placentae).

Coagulation / Fibrinolysis System

T01001-96-S SA,iQ,MX,A	Cardio Trombophilia Panel RT-PCR kit for detection of 8 mutations (F2 20210 G>A, F5 1691 G>A (Arg506Gln), F710976 G>A (Arg353Gln), F13 G>T (Val34Leu), FGB -455 G>A, ITGA2 807 C>T (Phe224Phe), ITGB3 1565 T>C(Leu33Pro), SERPINE1(PAI1) -675 5G>4G) in ready to use 12 x 8 strip format	R	96
T01002-96-S SA,iQ,MX,A	Folates Methabolism Panel RT-PCR kit for detection of 4 mutations: MTHFR 677 C>T (Ala222Val), MTHFR 1298A>C (Glu429Ala), MTR 2756 A>G (Asp3919Gly), MTRR 66 A>G (Ile22Met) in ready to use 12 x 8 strip format	R	96
T01101-50-T SA,RG,iQ,MX,A,B	FV (G1691A) Leiden SNP-Screen RT-PCR test for detection of Leiden mutation (Arg506Gln; rs6025) - ready to use 0,2 ml tube format	R	60
T01101-96-S SA,iQ,MX,A,B	FV (G1691A) Leiden SNP-Screen RT-PCR test for detection of Leiden mutation (Arg506Gln; rs6025) - ready to use 12 x 8 strip format	R	96
T01102-50-T SA,RG,iQ,MX,A,B	FII Protrombin (G20210A) SNP-Screen RT-PCR test for detection of prothrombin F2 gene mutation (rs1799963) - ready to use 0,2 ml tube format	R	60
T01102-96-S SA,RG,iQ,MX,A,B	FII Protrombin (G20210A) SNP-Screen RT-PCR test for detection of prothrombin F2 gene mutation (rs1799963) - ready to use 12 x 8 strip format	R	96
T01105-50-T SA,RG,iQ,MX,A,B	FVII (G1238A) SNP-Screen RT-PCR test for detection of FVII gene mutation (Arg353Gln; rs6046) - ready to use 0,2 ml tube format	R	60
T01105-96-S SA,RG,iQ,MX,A,B	FVII (G1238A) SNP-Screen RT-PCR test for detection of FVII gene mutation (Arg353Gln; rs6046) - ready to use 12 x 8 strip format	R	96
T01103-50-T SA,RG,iQ,MX,A,B	MTHFR (C677T) SNP-Screen RT-PCR test for detection of MTHFR gene mutation (Ala222Val; rs1801133) - ready to use 0,2 ml tube format	R	60
T01103-96-S SA,iQ,MX,A,B	MTHFR (C677T) SNP-Screen RT-PCR test for detection of MTHFR gene mutation (Ala222Val; rs1801133) - ready to use 12 x 8 strip format	R	96

Coagulation / Fibrinolysis System

T01273-50-T SA, RG, IQ, MX, A, B	MTHFR (A1298C) SNP-Screen RT-PCR test for detection of MTHFR gene mutation (Glu429Ala; Rs1801131) - ready to use 0,2 ml tube format	R	60
T01273-96-S SA, RG, IQ, MX, A, B	MTHFR (A1298C) SNP-Screen RT-PCR test for detection of MTHFR gene mutation (Glu429Ala; Rs1801131) - ready to use 12 x 8 strip format	R	96
T01124-50-T SA, RG, IQ, MX, A, B	MTRR (A 66 G) SNP-Screen RT-PCR test for detection of methioninesynthase gene mutation (Ile22Met; rs1801394) - ready to use 0,2 ml tube format	R	60
T01124-96-S SA, RG, IQ, MX, A, B	MTRR (A 66 G) SNP-Screen RT-PCR test for detection of methioninesynthase gene mutation (Ile22Met; rs1801394) - ready to use 12 x 8 strip format	R	96
T01143-50-T SA, RG, IQ, MX, A, B	MTR (A2756G) SNP-Screen RT-PCR test for detection of methionine synthase gene mutation (Asp919Gly; rs1805087) - ready to use 0,2 ml tube format	R	60
T01143-96-S SA, RG, IQ, MX, A, B	MTR (A2756G) SNP-Screen RT-PCR test for detection of methionine synthase gene mutation (Asp919Gly; rs1805087) - ready to use 12 x 8 strip format	R	96
T01120-50-T SA, RG, IQ, MX, A, B	PAI SERPINE (-675 5G/4G) SNP-Screen RT-PCR test for detection of insertion/deletion polymorphism of SERPINE1 or plasminogen activator inhibitor type 1 gene (rs1799768) - ready to use 0,2 ml tube format	R	60
T01120-96-S SA, RG, IQ, MX, A, B	PAI SERPINE (-675 5G/4G) SNP-Screen RT-PCR test for detection of insertion/deletion polymorphism of SERPINE1 or plasminogen activator inhibitor type 1 gene (rs1799768) - ready to use 12 x 8 strip format	R	96
T01107-50-T SA, RG, IQ, MX, A, B	FGB (G-455A) SNP-Screen RT-PCR test for detection of fibrinogen beta gene (rs1800790) - ready to use 0,2 ml tube format	R	60
T01107-96-S SA, RG, IQ, MX, A, B	FGB (G-455A) SNP-Screen RT-PCR test for detection of fibrinogen beta gene (rs1800790) - ready to use 12 x 8 strip format	R	96
T01356-50-T SA, RG, IQ, MX, A, B	FXII (C -4T) SNP-Screen RT-PCR test for detection of factor XII gene mutation (rs1801020) - ready to use 0,2 ml tube format	R	60
T01356-96-S SA, RG, IQ, MX, A, B	FXII (C -4T) SNP-Screen RT-PCR test for detection of factor XII gene mutation (rs1801020) - ready to use 12 x 8 strip format	R	96
T01355-50-T SA, RG, IQ, MX, A, B	FXIII (G103T) SNP-Screen RT-PCR test for detection of factor XIII A1 gene mutation (Val35Leu; rs5985) - ready to use 0,2 ml tube format	R	60
T01355-96-S SA, IQ, MX, A, B	FXIII (G103T) SNP-Screen RT-PCR test for detection of factor XIII A1 gene mutation (Val35Leu; rs5985) - ready to use 12 x 8 strip format	R	96

Coagulation / Fibrinolysis System

T01155-50-T SA, RG, IQ, MX, A, B	ITGA2 (C807T) SNP-Screen RT-PCR test for detection of Integrin alpha2 gene mutation (Phe224Phe; rs1126643) - ready to use 0,2 ml tube format	R	60
T01155-96-S SA, IQ, MX, A, B	ITGA2 (C807T) SNP-Screen RT-PCR test for detection of Integrin alpha2 gene mutation (Phe224Phe; rs1126643) - ready to use 12 x 8 strip format	R	96
T01106-50-T SA, RG, IQ, MX, A, B	ITGB3 (T176C) SNP-Screen RT-PCR test for detection of Integrin beta-3 gene mutation (Leu33Pro; rs5918) - ready to use 0,2 ml tube format	R	60
T01106-96-S SA, IQ, MX, A, B	ITGB3 (T176C) SNP-Screen RT-PCR test for detection of Integrin beta-3 gene mutation (Leu33Pro; rs5918) - ready to use 12 x 8 strip format	R	96
T01179-50-T SA, RG, IQ, MX, A, B	GPIBA1b (C482T) SNP-Screen RT-PCR test for detection of platelet glycoprotein Iba gene mutation (Thr145 Met; rs6065) - ready to use 0,2 ml tube format	R	60
T01179-96-S SA, IQ, MX, A, B	GPIBA1b (C482T) SNP-Screen RT-PCR test for detection of platelet glycoprotein Iba gene mutation (Thr145 Met; rs6065) - ready to use 12 x 8 strip format	R	96
T01354-50-T SA, RG, IQ, MX, A, B	GPIBA1b (T -5C) SNP-Screen RT-PCR test for detection of platelet glycoprotein Iba gene mutation (Kozak sequence polymorphism; rs2243093) - ready to use 0,2 ml tube format	R	60
T01354-96-S SA, IQ, MX, A, B	GPIBA1b (T -5C) SNP-Screen RT-PCR test for detection of platelet glycoprotein Iba gene mutation (Kozak sequence polymorphism; rs2243093) - ready to use 12 x 8 strip format	R	96
T01357-50-T SA, RG, IQ, MX, A, B	SELPLG (G186A) SNP-Screen RT-PCR test for detection selectin P ligand gene mutation (Met62Ile; rs2228315) - ready to use 0,2 ml tube format	R	60
T01357-96-S SA, IQ, MX, A, B	SELPLG (G186A) SNP-Screen RT-PCR test for detection selectin P ligand gene mutation (Met62Ile; rs2228315) - ready to use 12 x 8 strip format	R	96

Pharmacogenetics

Pharmacogenetics is the study of genetic differences in metabolic pathways which can affect individual responses to drugs, both in terms of therapeutic effect and adverse effects.

Warfarin (also known by the brand names Coumadin, Jantoven, Marevan, Uniwarfin) is an anticoagulant normally used in the prevention of thrombosis and thromboembolism, the formation of blood clots in the blood vessels and their migration elsewhere in the body, respectively. Polymorphisms in two genes, VKORC1 and CYP2C9, can affect the sensitivity of an individual patient to warfarin.

Warfarin Sensitivity

T01104-50-T SA, RG, IQ, MX, A, B	CYP2C9*2 (C430T) SNP-Screen RT-PCR test for detection of CYP2C9 gene mutation (Arg144Cys; rs1799853) - ready to use 0,2 ml tube format	R	60
T01104-96-S SA, IQ, MX, A, B	CYP2C9*2 (C430T) SNP-Screen RT-PCR test for detection of CYP2C9 gene mutation (Arg144Cys; rs1799853) - ready to use 12 x 8 strip format	R	96
T01111-50-T SA, RG, IQ, MX, A, B	CYP2C9*3 (A1075C) SNP-Screen RT-PCR test for detection of CYP2C9 gene mutation (Ile359Leu; rs1057910) - ready to use 0,2 ml tube format	R	60
T01111-96-S SA, IQ, MX, A, B	CYP2C9*3 (A1075C) SNP-Screen RT-PCR test for detection of CYP2C9 gene mutation (Ile359Leu; rs1057910) - ready to use 12 x 8 strip format	R	96
T01144-50-T SA, RG, IQ, MX, A, B	VKORC1 (C1173T) SNP-Screen RT-PCR test for detection Vitamin K epoxide reductase complex subunit 1 gene mutation, Warfarin sensitivity (rs9934438) - ready to use 0,2 ml tube format	R	60
T01144-96-S SA, IQ, MX, A, B	VKORC1 (C1173T) SNP-Screen RT-PCR test for detection Vitamin K epoxide reductase complex subunit 1 gene mutation, Warfarin sensitivity (rs9934438) - ready to use 12 x 8 strip format	R	96
T01145-50-T SA, RG, IQ, MX, A, B	VKORC1 (G3730A) SNP-Screen RT-PCR test for detection Vitamin K epoxide reductase complex subunit 1 gene mutation, Warfarin resistance (rs7294) - ready to use 0,2 ml tube format	R	60
T01145-96-S SA, IQ, MX, A, B	VKORC1 (G3730A) SNP-Screen RT-PCR test for detection Vitamin K epoxide reductase complex subunit 1 gene mutation, Warfarin resistance (rs7294) - ready to use 12 x 8 strip format	R	96




Fig. Warfarin mechanism of action

Clopidogrel (Plavix)

The anti-platelet agent clopidogrel bisulfate (sold under the trade name Plavix or Clopidogrel) is a widely prescribed medication for the prevention of blood clots in patients with acute coronary syndrome. This drug requires activation by CYP2C19, therefore individual carriers of alleles with reduced activity like CYP2C19*2 and CYP2C19*3 are likely at risk of therapeutic failure.

T01323-50-T SA, RG, IQ, MX, A, B	CYP2C19*2 (G681A) SNP-Screen RT-PCR test for detection of CYP2C19 gene mutation (rs4244285) - ready to use 0,2 ml tube format	R	60
T01323-96-S SA, IQ, MX, A, B	CYP2C19*2 (G681A) SNP-Screen RT-PCR test for detection of CYP2C19 gene mutation (rs4244285) - ready to use 12 x 8 strip format	R	96
T01324-50-T SA, RG, IQ, MX, A, B	CYP2C19*3 (G636A) SNP-Screen RT-PCR test for detection of CYP2C19 gene mutation (rs4986893) - ready to use 0,2 ml tube format	R	60
T01324-96-S SA, IQ, MX, A, B	CYP2C19*3 (G636A) SNP-Screen RT-PCR test for detection of CYP2C19 gene mutation (rs4986893) - ready to use 12 x 8 strip format	R	96

Tacrolimus (FK-506 or Fujimycin)

Tacrolimus, also known as FK-506 or Fujimycin, is the generic name for a calcineurin inhibitor drug. It is an immunosuppressive agent for treating autoimmune disease, including myasthenia gravis and rheumatoid arthritis, as well as for preventing allograft rejection in organ transplantation. Tacrolimus is primarily metabolized by CYP3A5. CYP3A5*3 is a nonfunctional variant, so heterozygous and especially homozygous carriers of this allele tend not to break down tacrolimus as much as normal allele, leading to higher blood concentrations of tacrolimus in these (CYP3A5*3) individuals.

T01331-50-T SA, RG, IQ, MX, A, B	CYP3A5*3 (G6986A) SNP-Screen RT-PCR test for detection of CYP3A5 gene mutation (rs776746) - ready to use 0,2 ml tube format	R	60
T01331-96-S SA, IQ, MX, A, B	CYP3A5*3 (G6986A) SNP-Screen RT-PCR test for detection of CYP3A5 gene mutation (rs776746) - ready to use 12 x 8 strip format	R	96

Statins

The SLCO1B1 gene encodes for the organic anion transporting polypeptide 1B1 (OATP1B1), an influx transporter produced in the liver that mediates the hepatic uptake and metabolism of statins.

Inherited variations in the SLCO1B1 gene known as SNPs (single nucleotide polymorphisms) affect the function of this transporter. The presence of this variant, especially in homozygotes, results in significantly decreased ability to take up statins, less effectiveness of the statin in lowering LDL-C, higher blood levels after dosing, and an increased risk of myopathy. Studies show that people who have particular inherited variations on the SLCO1B1 gene are four- to 17-times more likely to suffer myopathy as a side effect.

T01303-50-T SA, RG, IQ, MX, A, B	SLCO1B1 (T37041C) SNP-Screen RT-PCR test for detection of SLCO1B1 gene mutation (Val174Ala; rs4149056) - ready to use 0,2 ml tube format	R	60
T01303-96-S SA, IQ, MX, A, B	SLCO1B1 (T37041C) SNP-Screen RT-PCR test for detection of SLCO1B1 gene mutation (Val174Ala; rs4149056) - ready to use 12 x 8 strip format	R	96

Diabetes and obesity

Obesity is a medical condition in which excess body fat has accumulated to the extent that it may have a negative effect on health, leading to reduced life expectancy and/or increased health problems. Obesity increases the likelihood of various diseases, particularly heart disease, type 2 diabetes, obstructive sleep apnea, certain types of cancer, and osteoarthritis. Obesity is most commonly caused by a combination of excessive food energy intake, lack of physical activity, and genetic susceptibility.

Diabetes and obesity

T01335-50-T SA, RG, IQ, MX, A, B	PPARG2 (C34G) SNP-Screen Real Time PCR test for detection of Peroxisome Proliferator-Activated Receptor-Gamma-2 gene mutation (Pro12Ala, rs1801282) - ready to use 0,2 ml tube format	R	60
T01335-96-S SA, IQ, MX, A, B	PPARG2 (C34G) SNP-Screen Real Time PCR test for detection of Peroxisome Proliferator-Activated Receptor-Gamma-2 gene mutation (Pro12Ala, rs1801282) - ready to use 12 x 8 strip format	R	96
T01358-50-T SA, RG, IQ, MX, A, B	ADRB2 (C5318G) SNP-Screen RT-PCR test for detection adrenoceptor beta 2, surface gene (Gln27Glu, rs1042714) - ready to use 0,2 ml tube format	R	60
T01358-96-S SA, IQ, MX, A, B	ADRB2 (C5318G) SNP-Screen RT-PCR test for detection adrenoceptor beta 2, surface gene (Gln27Glu, rs1042714) - ready to use 12 x 8 strip format	R	96
T01359-50-T SA, RG, IQ, MX, A, B	ADRB2 (G46A) SNP-Screen RT-PCR test for detection adrenoceptor beta 2, surface gene (Arg16Gly, rs1042713) - ready to use 0,2 ml tube format	R	60
T01359-96-S SA, IQ, MX, A, B	ADRB2 (G46A) SNP-Screen RT-PCR test for detection adrenoceptor beta 2, surface gene (Arg16Gly, rs1042713) - ready to use 12 x 8 strip format	R	96
T01360-50-T SA, RG, IQ, MX, A, B	ADRB3 (T190C) SNP-Screen RT-PCR test for detection of β -adrenergic receptor genes mutation (Trp64Arg, rs4994) - ready to use 0,2 ml tube format	R	60
T01360-96-S SA, IQ, MX, A, B	ADRB3 (T190C) SNP-Screen RT-PCR test for detection of β -adrenergic receptor genes mutation (Trp64Arg, rs4994) - ready to use 12 x 8 strip format	R	96
T01361-50-T SA, RG, IQ, MX, A, B	FABP2 (A163G) SNP-Screen RT-PCR test for detection of Fatty Acid Binding Protein 2 gene (Ala54Thr, rs1799883) - ready to use 0,2 ml tube format	R	60
T01361-96-S SA, IQ, MX, A, B	FABP2 (A163G) SNP-Screen RT-PCR test for detection of Fatty Acid Binding Protein 2 gene (Ala54Thr, rs1799883) - ready to use 12 x 8 strip format	R	96
T01329-50-T SA, RG, IQ, MX, A, B	FTO (A23525T) SNP-Screen RT-PCR test for detection of FTO gene mutation (rs9939609) - ready to use 0,2 ml tube format	R	60
T01329-96-S SA, IQ, MX, A, B	FTO (A23525T) SNP-Screen RT-PCR test for detection of FTO gene mutation (rs9939609) - ready to use 12 x 8 strip format	R	96

Hepatitis C treatment prognosis

Chronic infection with hepatitis C virus (HCV) affects 170 million people worldwide and is the leading cause of cirrhosis in North America. Although the recommended treatment for chronic infection involves a 48-week course of peginterferon- α -2b (PegIFN- α -2b) or - α -2a (PegIFN- α -2a) combined with ribavirin (RBV), it is well known that many patients will not be cured by treatment.

Recent research has shown that genetic variation in the IL28B gene predicts both chronicity of HCV infection and sustained virological response (SVR) to antiviral standard therapy.

Single nucleotide polymorphisms (SNPs) near **interleukin-28B (IL-28B)** gene were shown to be highly associated with treatment response (SVR) in patients with chronic hepatitis C virus (HCV) infection.

Hepatitis C treatment prognosis

R05-100FRT SA, RG, iQ, SC, MX, A, B	IL28B rs17 / rs60 Real-TM RT-PCR test for detection of Interleukin mutations	R	100
T01349-50-T SA, RG, iQ, MX, A, B	IL28B (T>G) SNP-Screen RT-PCR test for detection of Interleukin mutation 1 (rs8099917) - ready to use 0,2 ml tube format	R	60
T01349-96-S SA, iQ, MX, A, B	IL28B (T>G) SNP-Screen RT-PCR test for detection of Interleukin mutation 1 (rs8099917) - ready to use 12 x 8 strip format	R	96
T01371-50-T SA, RG, iQ, MX, A, B	IL28B (C>T) SNP-Screen RT-PCR test for detection of Interleukin mutation 2 (rs12979860) - ready to use 0,2 ml tube format	R	60
T01371-96-S SA, iQ, MX, A, B	IL28B (C>T) SNP-Screen RT-PCR test for detection of Interleukin mutation 2 (rs12979860) - ready to use 12 x 8 strip format	R	96




Fig. 3D image of IL28B protein structure

Hypertension

Hypertension (HTN) or high blood pressure, sometimes called arterial hypertension, is a chronic medical condition in which the blood pressure in the arteries is elevated. Blood pressure is summarised by two measurements, systolic and diastolic, which depend on whether the heart muscle is contracting (systole) or relaxed between beats (diastole). As of 2000, nearly one billion people or ~26% of the adult population of the world had hypertension. It was common in both developed (333 million) and undeveloped (639 million) countries. In Europe hypertension occurs in about 30-45% of people as of 2013.

Genes most involved in Hypertension are the ones coding for angiotensin (**AGT**), Nitric oxide synthase 3 (**NOS3**) and Angiotensin Receptor (**AGTR**).

A genetic variant of the AGT gene leads to increased production of angiotensinogen. Therefore, carriers of this AGT variant have a higher risk for hypertension.

Nitric oxide is catalyzed by endothelial nitric oxide synthase (NOS), an enzyme with multiple genetic variants that might confer risk for hypertension. Single Nucleotide Polymorphism C786T has been associated with hypertension.

Polymorphism in the angiotensin II type 1 receptor (AGTR) gene is associated with the incidence of essential hypertension and increased coronary artery vasoconstriction.

Hypertension

T01118-50-T SA, RG, IQ, MX, A, B	AGT 1 (C521T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 1 (Thr207Met, rs4762) - ready to use 0,2 ml tube format	R	60
T01118-96-S SA, IQ, MX, A, B	AGT 1 (C521T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 1 (Thr207Met, rs4762) - ready to use 12 x 8 strip format	R	96
T01119-50-T SA, RG, IQ, MX, A, B	AGT 2 (C4072T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 2 (Met235Thr, rs699) - ready to use 0,2 ml tube format	R	60
T01119-96-S SA, IQ, MX, A, B	AGT 2 (C4072T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 2 (Met235Thr, rs699) - ready to use 12 x 8 strip format	R	96
T01182-50-T SA, RG, IQ, MX, A, B	NOS3 (C786T) SNP-Screen RT-PCR test for detection of Nitric oxide synthase 3 mutation (rs2070744) - ready to use 0,2 ml tube format	R	60
T01182-96-S SA, IQ, MX, A, B	NOS3 (C786T) SNP-Screen RT-PCR test for detection of Nitric oxide synthase 3 mutation (rs2070744) - ready to use 12 x 8 strip format	R	96
T01131-50-T SA, RG, IQ, MX, A, B	AGTR1 (A1166C) SNP-Screen RT-PCR test for detection of angiotensin II receptor, type 1 mutation (rs5186) - ready to use 0,2 ml tube format	R	60
T01131-96-S SA, IQ, MX, A, B	AGTR1 (A1166C) SNP-Screen RT-PCR test for detection of angiotensin II receptor, type 1 mutation (rs5186) - ready to use 12 x 8 strip format	R	96

Myocardial infarction / Ischemic stroke

Thrombotic complications, such as an acute myocardial infarction, ischemic stroke, pulmonary embolism, deep venous thrombosis are the major causes of morbidity and mortality in the world.

Myocardial infarction / Ischemic stroke

T01118-50-T SA, RG, IQ, MX, A, B	AGT 1 (C521T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 1 (Thr207Met, rs4762) - ready to use 0,2 ml tube format	R	60
T01118-96-S SA, IQ, MX, A, B	AGT 1 (C521T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 1 (Thr207Met, rs4762) - ready to use 12 x 8 strip format	R	96
T01119-50-T SA, RG, IQ, MX, A, B	AGT 2 (C4072T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 2 (Met235Thr, rs699) - ready to use 0,2 ml tube format	R	60
T01119-96-S SA, IQ, MX, A, B	AGT 2 (C4072T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 2 (Met235Thr, rs699) - ready to use 12 x 8 strip format	R	96
T01182-50-T SA, RG, IQ, MX, A, B	NOS3 (C786T) SNP-Screen RT-PCR test for detection of Nitric oxide synthase 3 mutation (rs2070744) - ready to use 0,2 ml tube format	R	60
T01182-96-S SA, IQ, MX, A, B	NOS3 (C786T) SNP-Screen RT-PCR test for detection of Nitric oxide synthase 3 mutation (rs2070744) - ready to use 12 x 8 strip format	R	96
T01148-50-T SA, RG, IQ, MX, A, B	APOE Leu48Pro SNP-Screen RT-PCR test for detection of apolipoprotein E mutation (T3100C, rs769452) - ready to use 0,2 ml tube format	R	60
T01148-96-S SA, IQ, MX, A, B	APOE Leu48Pro SNP-Screen RT-PCR test for detection of apolipoprotein E mutation (T3100C, rs769452) - ready to use 12 x 8 strip format	R	96
T01179-50-T SA, RG, IQ, MX, A, B	GP1BA1b (C482T) SNP-Screen RT-PCR test for detection of Platelet glycoprotein Iba mutation (Thr161Met, rs6065) - ready to use 0,2 ml tube format	R	60
T01179-96-S SA, IQ, MX, A, B	GP1BA1b (C482T) SNP-Screen RT-PCR test for detection of Platelet glycoprotein Iba mutation (Thr161Met, rs6065) - ready to use 12 x 8 strip format	R	96
T01354-50-T SA, RG, IQ, MX, A, B	GP1BA1b (T -5C) SNP-Screen RT-PCR test for detection of Platelet glycoprotein Iba mutation (rs2243093) - ready to use 0,2 ml tube format	R	60
T01354-96-S SA, IQ, MX, A, B	GP1BA1b (T -5C) SNP-Screen RT-PCR test for detection of Platelet glycoprotein Iba mutation (rs2243093) - ready to use 12 x 8 strip format	R	96
T01155-50-T SA, RG, IQ, MX, A, B	ITGA2 (C807T) SNP-Screen RT-PCR test for detection of integrin, alpha 2 mutation (Phe253Phe, rs1126643) - ready to use 0,2 ml tube format	R	60
T01155-96-S SA, IQ, MX, A, B	ITGA2 (C807T) SNP-Screen RT-PCR test for detection of integrin, alpha 2 mutation (Phe253Phe, rs1126643) - ready to use 12 x 8 strip format	R	96

Woman's Health / Miscarriage

Hypercoagulability in pregnancy, particularly due to inheritable thrombophilia, can lead to placental vascular thrombosis. This can in turn lead to complications like early-onset hypertensive disorders of pregnancy, pre-eclampsia and small for gestational age infants (SGA). Among other causes of hypercoagulability, Antiphospholipid syndrome has been associated with adverse pregnancy outcomes including recurrent miscarriage. Deep vein thrombosis has an incidence of one in 1,000 to 2,000 pregnancies in the United States, and is the second most common cause of maternal death in developed countries after bleeding.

Woman's Health / Miscarriage

T01101-50-T SA, RG, IQ, MX, A, B	FV (G1691A) Leiden SNP-Screen RT-PCR test for detection of Leiden mutation (Arg506Gln; rs6025) - ready to use 0,2 ml tube format	R	60
T01101-96-S SA, IQ, MX, A, B	FV (G1691A) Leiden SNP-Screen RT-PCR test for detection of Leiden mutation (Arg506Gln; rs6025) - ready to use 12 x 8 strip format	R	96
T01102-50-T SA, RG, IQ, MX, A, B	FII Protrombin (G20210A) SNP-Screen RT-PCR test for detection of prothrombin F2 gene mutation (rs1799963) - ready to use 0,2 ml tube format	R	60
T01102-96-S SA, RG, IQ, MX, A, B	FII Protrombin (G20210A) SNP-Screen RT-PCR test for detection of prothrombin F2 gene mutation (rs1799963) - ready to use 12 x 8 strip format	R	96
T01105-50-T SA, RG, IQ, MX, A, B	FVII (G1238A) SNP-Screen RT-PCR test for detection of FVII gene mutation (Arg353Gln; rs6046) - ready to use 0,2 ml tube format	R	60
T01105-96-S SA, RG, IQ, MX, A, B	FVII (G1238A) SNP-Screen RT-PCR test for detection of FVII gene mutation (Arg353Gln; rs6046) - ready to use 12 x 8 strip format	R	96
T01103-50-T SA, RG, IQ, MX, A, B	MTHFR (C677T) SNP-Screen RT-PCR test for detection of MTHFR gene mutation (Ala222Val; rs1801133) - ready to use 0,2 ml tube format	R	60
T01103-96-S SA, IQ, MX, A, B	MTHFR (C677T) SNP-Screen RT-PCR test for detection of MTHFR gene mutation (Ala222Val; rs1801133) - ready to use 12 x 8 strip format	R	96



Woman's Health / Miscarriage

T01273-50-T SA, RG, iQ, MX, A, B	MTHFR (A1298C) SNP-Screen RT-PCR test for detection of MTHFR gene mutation (Glu429Ala; Rs1801131) - ready to use 0,2 ml tube format	R	60
T01273-96-S SA, RG, iQ, MX, A, B	MTHFR (A1298C) SNP-Screen RT-PCR test for detection of MTHFR gene mutation (Glu429Ala; Rs1801131) - ready to use 12 x 8 strip format	R	96
T01124-50-T SA, RG, iQ, MX, A, B	MTRR (A 66 G) SNP-Screen RT-PCR test for detection of methioninesynthase gene mutation (Ile22Met; rs1801394) - ready to use 0,2 ml tube format	R	60
T01124-96-S SA, RG, iQ, MX, A, B	MTRR (A 66 G) SNP-Screen RT-PCR test for detection of methioninesynthase gene mutation (Ile22Met; rs1801394) - ready to use 12 x 8 strip format	R	96
T01143-50-T SA, RG, iQ, MX, A, B	MTR (A2756G) SNP-Screen RT-PCR test for detection of methionine synthase gene mutation (Asp919Gly; rs1805087) - ready to use 0,2 ml tube format	R	60
T01143-96-S SA, RG, iQ, MX, A, B	MTR (A2756G) SNP-Screen RT-PCR test for detection of methionine synthase gene mutation (Asp919Gly; rs1805087) - ready to use 12 x 8 strip format	R	96
T01120-50-T SA, RG, iQ, MX, A, B	PAI SERPINE (-675 5G/4G) SNP-Screen RT-PCR test for detection of insertion/deletion polymorphism of SERPINE1 or plasminogen activator inhibitor type 1 gene (rs1799768) - ready to use 0,2 ml tube format	R	60
T01120-96-S SA, RG, iQ, MX, A, B	PAI SERPINE (-675 5G/4G) SNP-Screen RT-PCR test for detection of insertion/deletion polymorphism of SERPINE1 or plasminogen activator inhibitor type 1 gene (rs1799768) - ready to use 12 x 8 strip format	R	96
T01107-50-T SA, RG, iQ, MX, A, B	FGB (G-455A) SNP-Screen RT-PCR test for detection of fibrinogen beta gene (rs1800790) - ready to use 0,2 ml tube format	R	60
T01107-96-S SA, RG, iQ, MX, A, B	FGB (G-455A) SNP-Screen RT-PCR test for detection of fibrinogen beta gene (rs1800790) - ready to use 12 x 8 strip format	R	96
T01118-50-T SA, RG, iQ, MX, A, B	AGT 1 (C521T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 1 (Thr207Met, rs4762) - ready to use 0,2 ml tube format	R	60
T01118-96-S SA, iQ, MX, A, B	AGT 1 (C521T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 1 (Thr207Met, rs4762) - ready to use 12 x 8 strip format	R	96
T01119-50-T SA, RG, iQ, MX, A, B	AGT 2 (C4072T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 2 (Met235Thr, rs699) - ready to use 0,2 ml tube format	R	60
T01119-96-S SA, iQ, MX, A, B	AGT 2 (C4072T) SNP-Screen RT-PCR test for detection of Angiotensinogen mutation 2 (Met235Thr, rs699) - ready to use 12 x 8 strip format	R	96

Immunity

TNF polymorphism influences immune function and regulation and moreover has been reported to be associated with Non-Hodgkin Lymphoma (NHL) risk. IL17A polymorphism is associated with a significant increase risk for specific type of cancer, especially in gastric cancer.

Immunity

T01177-50-T SA, RG, IQ, MX, A, B	TNF (G-308A) SNP-Screen RT PCR test for detection of Tumor necrosis factor mutation (rs1800629) - ready to use 0,2 ml tube format	R	60
T01177-96-S SA, IQ, MX, A, B	TNF (G-308A) SNP-Screen RT PCR test for detection of Tumor necrosis factor mutation (rs1800629) - ready to use 12 x 8 strip format	R	96
T01171-50-T SA, RG, IQ, MX, A, B	IL17A (G-197A) SNP-Screen RT PCR test for detection of Interleukin-17A gene mutation (rs2275913) - ready to use 0,2 ml tube format	R	60
T01171-96-S SA, IQ, MX, A, B	IL17A (G-197A) SNP-Screen RT PCR test for detection of Interleukin-17A gene mutation (rs2275913) - ready to use 12 x 8 strip format	R	96

Lipid Metabolism

Lipid metabolism refers to the processes that involve the intercourse and degradation of lipids. Certain abnormalities in these enzymes can lead to the buildup of specific fatty substances that normally would have been broken down by the enzymes. Over time, accumulations of these substances can be harmful to many organs of the body.

Lipid Metabolism

T01148-50-T SA, RG, IQ, MX, A, B	APOE Leu48Pro SNP-Screen RT PCR test for detection of Apolipoprotein E gene mutation (T3100C, rs769452) - ready to use 0,2 ml tube format	R	60
T01148-96-S SA, IQ, MX, A, B	APOE Leu48Pro SNP-Screen RT PCR test for detection of Apolipoprotein E gene mutation (T3100C, rs769452) - ready to use 12 x 8 strip format	R	96
T01149-50-T SA, RG, IQ, MX, A, B	LPL (Ser447Ter) SNP-Screen RT PCR test for detection of lipoprotein lipase gene mutation (1421C>G, rs328) - ready to use 0,2 ml tube format	R	60
T01149-96-S SA, IQ, MX, A, B	LPL (Ser447Ter) SNP-Screen RT-PCR test for detection of CYP3A5 gene mutation (rs776746) - ready to use 12 x 8 strip format	R	96
T01125-50-T SA, RG, IQ, MX, A, B	PON1 (Gln192Arg) SNP-Screen RT PCR test for detection of paraoxonase 1 gene mutation (575A>G, rs662) - ready to use 0,2 ml tube format	R	60
T01125-96-S SA, IQ, MX, A, B	PON1 (Gln192Arg) SNP-Screen RT PCR test for detection of paraoxonase 1 gene mutation (575A>G, rs662) - ready to use 12 x 8 strip format	R	96
T01161-50-T SA, RG, IQ, MX, A, B	LIPC (-250 G>A) SNP-Screen RT PCR test for detection of Hepatic lipase gene mutation (rs2070895) - ready to use 0,2 ml tube format	R	60
T01161-96-S SA, IQ, MX, A, B	LIPC (-250 G>A) SNP-Screen RT PCR test for detection of Hepatic lipase gene mutation (rs2070895) - ready to use 12 x 8 strip format	R	96

Men's Health

The Y chromosome accumulates male-related genes including sex-determining region of Y-chromosome (SRY) and several spermatogenesis-related genes. The long arm contains azoospermia factor (AZF) region (including sub-regions AZFa, AZFb and AZFc). Microdeletions in this region are responsible for azoospermia and oligospermia and result in the male infertility.

Men's Health

01200-50 SA,RG,IQ,MX,A,B	AZF System Y-chromosome RT-PCR test for detection of the microdeletions in AZF regions: AZFa (sY84, sY86), AZFb (sY127, sY134), AZFc (sY254, sY255) of the human Y chromosome	R	60
-----------------------------	--	---	----



HLA

Human Leukocyte Antigen (HLA) B27 (subtypes B*2701-2759) is a class I surface antigen encoded by the B locus in the major histocompatibility complex (MHC) on chromosome 6 and presents antigenic peptides (derived from self and non-self antigens) to T cells. HLA-B27 is strongly associated with ankylosing spondylitis (AS), and other associated inflammatory diseases referred to collectively as "spondyloarthritis".

HLA

R116-50FRT SA,RG,IQ,SC,MX,A,B	HLA B27 Real-TM Real Time PCR kit	R	48
----------------------------------	---	---	----





Fig. 3D image of HLA B27 protein structure

Oncological diseases

Chronic myelogenous leukaemia (CML) results from the neoplastic transformation of a haematopoietic stem cell. The hallmark genetic abnormality of CML is a t(9;22)(q34;q11) translocation, which was first discovered as an abnormal, small chromosome, named the 'Philadelphia chromosome'. This translocation generates the BCR - ABL fusion gene.

The discovery that BCR-ABL is required for the pathogenesis of CML, and that the tyrosine-kinase activity of ABL is essential for BCR-ABL-mediated transformation, made the ABL kinase an attractive target for therapeutic intervention. Imatinib mesylate (Glivec, previously known as STI571 and CGP 57148) — a potent inhibitor of the tyrosine kinases ABL, ARG, platelet-derived growth factor receptor and KIT — has been shown to selectively induce apoptosis of BCR-ABL+ cells, and is remarkably successful in treating patients with CML.




CML Kits


TR-O1 SA, RG, IQ, MX, SC, A, B	Mbcr-abl FRT Real-TM Complete Real Time Test with Ribo-Zol-D extraction kit	R	€	100
R-O1 SA, RG, IQ, MX, SC, A, B	Mbcr-abl FRT Real-TM Real Time Amplification kit	R	€	100

BRCA mutations are involved in development **breast and ovarian cancer**. BRCA1 and BRCA2 are tumor suppressor genes that are inactivated during neoplastic development in breast cancer. Germline mutations of the two genes are transmitted in the autosomal dominant way and predispose carriers to the development of ovarian and/or breast cancers. Mutations in BRCA1 are present in approximately one-half of the early-onset breast cancer families and 80% of the early-onset breast and ovarian cancer families, whereas BRCA2 mutations are believed to account for a comparable percentage of inherited breast cancer cases. Women with germline mutations in BRCA1 have a lifetime risk of 85% and up to 50% for breast and ovarian cancers, respectively. Approximately 50% to 65% of women born with a deleterious mutation in BRCA1 will develop breast cancer by age 70, and 35% to 46% will develop ovarian cancer by age 70. Approximately 40% to 57% of women with a deleterious mutation in BRCA2 will develop breast cancer by age 70, and 13% to 23% will develop ovarian cancer by age 70. Women with a breast cancer associated with a BRCA mutation have up to a 40% probability of developing a new primary breast cancer within 10 years following initial diagnosis if they did not receive tamoxifen treatment or have an oophorectomy. The woman's ten-year risk for ovarian cancer is also increased by 6-12% under these conditions.

BRCA mutations frequency in breast cancer (caucasian population)



BRCA mutations frequency in ovarian cancer (caucasian population)



BRCA Kits

R-27/P-48FRT SA, MX, IQ, B, A	Oncogenetics BRCA Panel Real-TM Real Time Amplification kit	R		48
----------------------------------	---	---	--	----

Cancer therapy

With the development of **molecular cancer biology**, targeted cancer therapy got much attention and selecting appropriate patients for targeted therapies became the key to treatment success. The epidermal growth factor receptor (**EGFR**) has become an important target for cancer treatment. Further studies have shown that the **mutational status** of **KRAS** and **BRAF** genes have a close association with the efficacy of EGFR targeted therapy. Patients carried KRAS, or BRAF mutations in their tumors had significant lower response rate to EGFR targeted therapy. The mutational status of **PIK3CA** gene has a close association with the efficacy of EGFR or HER2 targeted therapy.

Cancer therapy

BT001 SA,iQ,MX,A,B	Human KRAS Gene Real-TM NEW Real Time PCR test for detection of KRAS gene mutations: Exon 12, 34G>T mutation (G12C); 34G>A mutation (G12S); 34G>C mutation (G12R); 35G>T mutation (G12V); 35G>A mutation (G12D); 35G>C mutation (G12A); Exon 13, 38G>A mutation (G13D); ready to use 12 x 8 strip format	R	12
BT002 SA,iQ,MX,A,B	Human EGFR Gene Real-TM NEW Real Time PCR test for detection of more than 30 mutations in EGFR gene on Exons 18, 19, 20, 21 - ready to use 12 x 8 strip format	R	12
TBT004 SA,iQ,MX,A,B	Human BRAF Gene Real-TM NEW Real Time PCR test for detection of BRAF gene mutation GTG>GAG mutation(V600E) - ready to use 12 x 8 strip format	R	24
BT005 SA,iQ,MX,A,B	Human CKIT Gene D816V Real-TM NEW Real Time PCR test for detection of CKIT gene D816V mutation	R	24
BT006 SA,iQ,MX,A,B	Human PIK3CA Gene Real-TM NEW Real Time PCR test for detection of 9 PIK3CA gene mutations	R	24
BT007 SA,iQ,MX,A,B	Human KRAS Gene Codon 61 Real-TM NEW Real Time PCR test for detection of KRAS gene codon 61 mutation	R	12
BT008 SA,iQ,MX,A,B	Human NRAS Gene Real-TM NEW Real Time PCR test for detection of NRAS gene mutations	R	12
RH001 SA,iQ,MX,A,B	Human EML4-ALK Gene Fusion Real-TM NEW Real Time PCR test for detection of EML4-ALK fusion gene	R	12
RH002 SA,iQ,MX,A,B	Human ROS1 Gene Fusion Real-TM NEW Real Time PCR test for detection of ROS1 fusion gene	R	12
RH003 SA,iQ,MX,A,B	Human KIF5B-RET Gene Fusion Real-TM NEW Real Time PCR test for detection of KIF5B-RET fusion gene	R	12

Periodontitis

Periodontitis is a serious gum infection that damages the soft tissue and destroys the bone that supports the tooth. Periodontitis can cause tooth loss or worse, an increased risk of heart attack or stroke and other serious health problems. Bacterial plaque is believed to be the principal etiological factor in the onset and progression of periodontitis.

Porphyromonas gingivalis and *Tannerella forsythia* are strong markers of periodontitis in adults, and these species have been linked to the progression of the disease.

Periodontitis Kit

T01707-96-S SA,iQ,MX,A	PeriodontScreen Real-TM RT-PCR kit for detection of <i>Porphyromonas endodontalis</i> , <i>Porphyromonas gingivalis</i> , <i>Aggregatibacter actinomycetemcomitans</i> , <i>Treponema denticola</i> , <i>Fusobacterium nucleatum</i> , <i>Prevotella intermedia</i> , <i>Tannerella forsythia</i> - ready to use 12 x 8 strip format	R	96
---------------------------	---	---	----

Automatic Extraction System

SM003	SaMag Viral Nucleic Acids Extraction kit For use with SaMag-12 or SaMag-24 instruments; extraction of Viral nucleic acids (DNA/RNA) from human biological specimens such as serum, plasma, and other cell-free fluids	€€	48
SM016	SaMag Viral Nucleic Acids Large Volume Extraction kit NEW For use with SaMag-12 or SaMag-24 instruments; extraction of Viral nucleic acids (DNA/RNA) from human biological specimens such as plasma, serum or cell-free body fluids starting from up to 1000 µl of sample .	€€	48
SM006	SaMag Bacterial DNA Extraction kit For use with SaMag-12 or SaMag-24 instruments; extraction of genomic DNA from both Gram-positive and Gram-negative bacteria starting from various sample types	€€	48
SM007	SaMag STD DNA Extraction kit For use with SaMag-12 or SaMag-24 instruments; extraction of STD DNA from human biological specimens such as cervical, urethral, conjunctival swabs, urine sediment, prostatic liquid, seminal liquid	€€	48
SM008	SaMag TB DNA Extraction kit For use with SaMag-12 or SaMag-24 instruments; extraction of genomic DNA of Mycobacteria spp. (e.g. Mycobacterium tuberculosis) from clinical specimens or cultures	€€	48
SM004	SaMag Tissue DNA Extraction kit For use with SaMag-12 or SaMag-24 instruments; extraction of genomic DNA from a variety of tissues	€€	48
SM001	SaMag Blood DNA Extraction kit For use with SaMag-12 or SaMag-24 instruments; extraction of genomic DNA from whole blood, peripheral blood mononuclear cells or buffy coat	€€	48
SM009	SaMag FFPE DNA Extraction kit For use with SaMag-12 or SaMag-24 instruments; extraction of genomic DNA from FFPE samples	€€	48
SM015	SaMag Total RNA/DNA Extraction kit NEW For use with SaMag-12 or SaMag-24 instruments; extraction of total RNA/DNA from whole blood, blood cells, animal tissue, plant tissue, yeast or cultured cells samples	€€	48

Silica Sorbtion Method

K-1-1/A	DNA-Sorb-A DNA Extraction Kit from clinical material (smears, scrapes, urine)	€€	50
K-1-1/B	DNA-Sorb-B DNA extraction kit from whole blood, plasma, liquor, sputum, biopats, fecal extract, etc.	€€	50
K-1-6/50	DNA-Sorb-C DNA extraction kit from plants material, biopats, human tissues, food samples		50
K-2-1	Ribo-Sorb RNA/DNA purification kit from plasma, serum, liquor, tissue, feces	€€	50
K-2-1/100	Ribo-Sorb RNA/DNA purification kit from plasma, serum, liquor, tissue, feces, etc		100

Magnetic Beads RNA/DNA Extraction Method

K-2-16/1000	Magno-Virus For rapid magnetic beads purification of viral RNA and DNA (e.g. HCV, HIV, HBV, HAV, HDV, Enteroviruses, CMV) from cell free body fluids such as plasma or serum and for pools screening of human plasma comprised of equal aliquots of not more than 10 individual specimens (1 ml of pool composed from 10 samples).	€€	100
-------------	---	----	-----

Heat-based Thermic Method

K-2-17	Express-DNA-A Fast heat-based thermic DNA extraction from urogenital swabs, throat swabs, conjunctiva swabs, erosive and ulcerative elements of mucous membranes and skin, urine.	50
--------	---	----

Chemical-based Method

K-2-9	DNA/RNA Prep Guanidine/isopropanol RNA/DNA extraction kit from clinical materials: peripheral blood, plasma, cerebrospinal fluid, amniotic fluid, sputum, swabs, tissue, feces, etc	50
-------	---	----

Spin Column Purification Method

K-2/C	Ribo Virus CC For rapid purification of viral RNA and DNA (HCV, HBV, HAV, HIV, CMV, HSV, VZV, EBV, parvovirus B19, H5N1) from cell-free biological fluids	50
-------	---	----

K-1-1/E	Genomic column DNA Express Genomic DNA from whole blood, serum, plasma, buffy coat, platelets, body fluids	50
---------	--	----



Step	Description
1. Lyse	 600 µl RAV1 150 µl sample 70°C 5 min
2. Adjust DNA	600 µl Ethanol
3. Bind	 Load sample  1 min 8000 x g
4. Wash	 1st wash: 500 µl RAW 2nd wash: 600 µl RAV3 3rd wash: 200 µl RAV3  1st and 2nd 1 min 8000 x g  3rd 5 min 11000 x g
5. Elute	 50 µl Rnase-free H₂O or Buffer RE (70°C) 1-2 min  1 min 11000 x g

Table. RNA/DNA RiboVirus short protocol


SaCycler-96 Real Time PCR System **CE - IVD**

Features


- 5 fluorescence channels multiplexing for discrimination of up to five targets in a single reaction well
- Strong flexibility thanks to the 96-well format suited for standard PCR microplates, test tubes and strips
- Possibility to use tubes with different heights with regulation of the matrix block
- Optimal signal/noise ratio and absence of crosstalk ensured by the unique design of the optical track including a separate light source for each channel and a matrix CCD camera
- Light emitting diodes (LED) as a light source with a lifetime of about 100,000 hours that does not require maintenance or constant monitoring
- Wide dynamic range of detection using multiple exposure method, which leads the optimization of signal registration conditions to a whole new level, greatly simplifying or even eliminating the need for fluorescence settings
- Main applications are Real-Time quantitation, single nucleotide polymorphisms (SNPs) genotyping, melting curve and gene expression analysis
- Instrument comes with the most used Sacace protocols already installed in the software minimizing possibility of error during programming.



Quantitative analysis with the use of calibration standards, allows to determine the quantity of DNA target in the sample



Qualitative analysis allows to determine presence or absence of target DNA in a sample



Melting curves are applied to analysis for determination of polymorphisms of single nucleotide

Product code	Description	
SC-96I	SaCycler-96 Real Time PCR System IVD Real Time PCR Thermalcycler closed system with 5 fluorescence channels	CE
SC-96R	SaCycler-96 Real Time PCR System Real Time PCR Thermalcycler open system with 5 fluorescence channels	

SaMag - 12 / 24™ Automatic Nucleic Acids Extraction System

Features

- Fully automated high yield nucleic acid extraction with magnetic beads and PCR setup (for Sacace kits like HCV, HBV, HIV Dx)
- Reagent identification with on-board barcode scanners
- Pre-programmed protocols
- Flexible batch size from 1 up to 24 samples in parallel
- All required tubes, tips, plastics are inside the provided extraction kits
- Ready to use reagents
- Very fast extraction protocol (~ 40 minutes)
- Very simple operation (easy to install, operate, maintain)
- Isolation of very pure nucleic acids



SaMag - 12™ CE - IVD

SaMag-12™ system is a compact bench-top extractor for automated nucleic acid purification. It is able to process 1-12 samples in parallel, recommended for labs with small-medium sample throughputs



SaMag - 24™ CE - IVD

SaMag-24™ system doubles the sample capacity of the SaMag-12™ with 2 independently programmable modules for purification of up to 24 samples and maximum flexibility. It is recommended for labs with medium-high sample throughputs

Product code	Description
SM-12	SaMag-12 Automatic Nucleic Acid Extraction System for processing up to 12 samples per run
SM-24	SaMag-24 Automatic Nucleic Acid Extraction System with two independent modules for processing up to 24 samples per run

Terms and Conditions

VALIDITY:

This catalogue is valid from the 1st of November 2015 and invalidates all prior publications.

PRICE TERMS:

EXW Sacace company Como (Italy) in EUR and in accordance with Incoterms 2000.

ORDERS:

Orders may be placed by fax (+39 0314892926) or by e-mail (order@sacace.com) using the order form downloadable from www.sacace.com

DELIVERY:

The general delivery time for all our products is 3-4 weeks from receipt of order/date of SACACE's official Order Confirmation. The delivery of Sacace Biotechnologies Srl products will be by the best means in the judgement of Sacace Biotechnologies Srl unless otherwise specified.

NOTE: delivery time is understood to be the period of time that passes from receipt of prepayment in SACACE's bank account or receipt of irrevocable, documentary, confirmed Letter of Credit until shipment of goods from SACACE company.

COMPLAINTS:

All complaints must be submitted in writing to Sacace Biotechnologies Srl within 15 days from receipt of the goods.

USES:

The products sold by Sacace Biotechnologies Srl are laboratory reagents and are for research & diagnostics use only. The products must be stored as recommended in the individual product manuals supplied with each product. All products must be used by qualified personnel. Sacace Biotechnologies Srl will not be liable for any incidental, consequential or contingent damages involving the use of his products.



Sacace Biotechnologies Srl

via Scalabrini, 44. – 22100 –Como – Italy Tel +390314892927 Fax +390314892926 VAT: 01294510621
mail: info@sacace.com web: www.sacace.com



www.sacace.com

info@sacace.com

COMO, ITALY